

## Datasheet: MCA2526

**BATCH NUMBER 160588**

<b>Description:</b>	MOUSE ANTI HUMAN D-DIMER
<b>Specificity:</b>	D-dimer
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	DD5
<b>Isotype:</b>	IgG2b
<b>Quantity:</b>	0.2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting (1)	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1)This antibody recognises human D-dimer under non-reducing conditions.**

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein</b>	IgG concentration 1.0mg/ml

## Concentrations

---

**Immunogen** Homogenised fibrin clot.

---

## External Database Links

### UniProt:

[P02671](#) [Related reagents](#)

[P02675](#) [Related reagents](#)

[P02679](#) [Related reagents](#)

### Entrez Gene:

[2243](#) FGA [Related reagents](#)

[2244](#) FGB [Related reagents](#)

[2266](#) FGG [Related reagents](#)

---

**RRID** AB\_609621

---

**Fusion Partners** Spleen cells from immunised Balb/c mice were fused with cells of the Sp2/0 myeloma cell line.

---

**Specificity** **Mouse anti Human D-dimer antibody, clone DD5** recognises human D-dimer and shows low cross-reactivity with fibrinogen (D-monomer).

D-dimers are the principle degradation products of cross-linked fibrin, released into the bloodstream following the digestion of fibrin clots (thrombus) by the enzyme plasmin, and thus measurement of the level of D-dimer in plasma, can be used as an indication of thrombotic disorders. Conditions in which levels of D-dimer are elevated include deep vein thrombosis (DVT), venous thromboembolism (VTE), cardiovascular disease and disseminated intravascular coagulation (DIC), a possible acute condition arising from surgical operations, septic shock and liver disease.

---

## References

1. Rowbotham, B.J. *et al.* (1987) Measurement of crosslinked fibrin derivatives--use in the diagnosis of venous thrombosis. [Thromb Haemost. 57 \(1\): 59-61.](#)
  2. Kruskal, J.B. *et al.* (1992) Elevated fibrin-related and fibrinogen-related antigens in patients with liver disease. [Hepatology. 16 \(4\): 920-3.](#)
- 

## Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

**Guarantee** 12 months from date of despatch

---

## Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2526>  
10040

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M389435:210806'

**Printed on 18 Jan 2024**

---