

Datasheet: MCA2559A647

Description:	MOUSE ANTI HUMAN CD263:Alexa Fluor® 647
Specificity:	CD263
Other names:	DcR1, TRAIL-R3
Format:	ALEXA FLUOR® 647
Product Type:	Monoclonal Antibody
Clone:	2B8
Isotype:	IgG1
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Alexa Fluor® 647 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®647	650	665
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide (NaN ₃)		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.05mg/ml		

Immunogen	DcR1-Fc chimaeric protein.
External Database Links	<p>UniProt: O14798 Related reagents</p> <p>Entrez Gene: 8794 TNFRSF10C Related reagents</p>
Synonyms	DCR1, LIT, TRAILR3, TRID
RRID	AB_906152
Fusion Partners	Lymph node cells from immunized mice were fused with cells of the P3U1 myeloma cell line
Specificity	<p>Mouse anti Human CD263 antibody, clone 2B8 recognizes human CD263, a ~60 kDa cell surface protein, which is otherwise known as decoy receptor 1 (DcR1). CD263 is a member of the tumour necrosis factor (TNF) receptor superfamily that binds to TNF-related apoptosis inducing ligand (TRAIL). As the receptor does not possess cytoplasmic death domains it cannot induce apoptosis, but instead plays an important role in inhibiting TRAIL-induced apoptosis. CD263 is expressed in normal tissues including heart, skeletal tissues and peripheral blood neutrophils, but weakly expressed in tumour cells.</p> <p>The 2B8 antibody does not inhibit DcR1 (CD263) inhibitory function in TRAIL-induced apoptosis.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Acknowledgements	This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad

Health And Safety Information Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:Alexa Fluor® 647 \(MCA928A647\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M367178:200529'

Printed on 04 Jan 2021