

Datasheet: MCA2559PE

Description:	MOUSE ANTI HUMAN CD263:RPE
Specificity:	CD263
Other names:	DcR1, TRAIL-R3
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	2B8
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1.0ml distilled water.		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide (NaN ₃)		
Stabilisers	1%	Bovine Serum Albumin	
	5%	Sucrose	

Immunogen	DcR1-Fc chimaeric protein.
External Database Links	<p>UniProt: O14798 Related reagents</p> <p>Entrez Gene: 8794 TNFRSF10C Related reagents</p>
Synonyms	DCR1, LIT, TRAILR3, TRID
RRID	AB_961534
Fusion Partners	Lymph node cells from immunized mice were fused with cells of the P3U1 myeloma cell line
Specificity	<p>Mouse anti Human CD263 antibody, clone 2B8 recognizes human CD263, a ~60 kDa cell surface protein, which is otherwise known as decoy receptor 1 (DcR1). CD263 is a member of the tumour necrosis factor (TNF) receptor superfamily that binds to TNF-related apoptosis inducing ligand (TRAIL). As the receptor does not possess cytoplasmic death domains it cannot induce apoptosis, but instead plays an important role in inhibiting TRAIL-induced apoptosis. CD263 is expressed in normal tissues including heart, skeletal tissues and peripheral blood neutrophils, but weakly expressed in tumour cells.</p> <p>The 2B8 antibody does not inhibit DcR1 (CD263) inhibitory function in TRAIL-induced apoptosis.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20487 available at: 20487: https://www.bio-rad-antibodies.com/uploads/MSDS/20487.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M375517:210104'

Printed on 23 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)