

Datasheet: MCA2600

Description:	MOUSE ANTI PIG GRANULOCYTES
Specificity:	GRANULOCYTES
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	2B2
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

RRID AB_1102598

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/100 - 1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Pig

Product Form Purified IgG - liquid

Preparation Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Porcine bone marrow hematopoietic cells (BMHC)

Fusion Partners Spleen cells from immunized Balb/c mouse were fused with cells of the Sp2/0 mouse myeloma cell line

Specificity **Mouse anti Pig granulocytes antibody, clone 2B2** recognizes porcine basophils, neutrophils and eosinophils, acting as a reliable tool for their analysis and isolation, without contamination from other cells.

Analysis of porcine bone marrow and blood shows that the antigen recognized by clone 2B2 is selectively expressed on mature cells and that it is likely to be a labile conformational epitope. Use of clone 2B2 in conjunction with clone 6D10 product code ([MCA2599GA](#)), allows for the discrimination and characterization of different porcine granulocyte lineages and also their developmental stages: 6D10⁻2B2⁻ early myeloid precursors, 6D10⁺2B2⁻ immature neutrophils, 6D10⁺2B2⁺ mature neutrophils and 6D10⁻2B2⁺ mature eosinophils and basophils. Clone 2B2 has been shown as suitable for use on cytopins ([Perez et al. 2007](#)).

Flow Cytometry Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.

References

1. Pérez, C. *et al.* (2007) Phenotypic and functional characterization of porcine granulocyte developmental stages using two new markers. [Dev Comp Immunol. 31 \(3\): 296-306.](#)
2. Nguyen, D.N. *et al.* (2016) Delayed development of systemic immunity in preterm pigs as a model for preterm infants. [Sci Rep. 6: 36816.](#)

Further Reading

1. Piriou-guzylack, L. & Salmon, H. (2008) Membrane markers of the immune cells in swine: an update. [Vet Res. 39 \(6\): 54.](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®549](#),
[DyLight®649](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com
'M334687:181203'

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

Printed on 02 Jan 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)