

## Datasheet: MCA2611

**BATCH NUMBER 160271**

<b>Description:</b>	MOUSE ANTI HUMAN SC5b-9
<b>Specificity:</b>	SC5b-9
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	056B-75.2.3.10 (3R2/0)
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen	▪			
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Borate buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	Current, batch-specific concentration 1.10 mg/ml
<b>RRID</b>	AB_844532

<b>Specificity</b>	<p><b>Mouse anti Human SC5b-9 antibody, clone 056B-75.2.3.10 (3R2/0)</b> recognizes the SC5b-9 complex of around 330 kDa. Both the classical and alternative complement pathways result in the formation of the cytolysis inducing C5b-9 complex. This complex is composed of 190 kDa C5b which is bound to cells via 71 kDa C9. Sublytic assembly of C5b-9 on plasma membranes induces cell cycle activation and survival. The binding of C5b-9 to the 75 kDa S-protein (or vitronectin) in the fluid phase prevents C5b-9 from assembling on the plasma membrane, deactivating it and forming the SC5b-9 complex. SC5b-9 is stable in vitro and is therefore a reliable indicator of terminal complement pathway activation.</p> <p>Removal of Sodium Azide is recommended prior to use in functional assays.</p> <p>Mouse anti Human SC5b-9 antibody, clone 3R2/0 has been reported to work in western blotting applications.</p>
<b>Histology Positive Control Tissue</b>	Kidney from post streptococcal glomerulonephritis patients.
<b>Further Reading</b>	<ol style="list-style-type: none"> <li>1. Biesecker, G. (1990) The complement SC5b-9 complex mediates cell adhesion through a vitronectin receptor. <a href="#">J Immunol. 145 (1): 209-14.</a></li> <li>2. Greenstein, J.D. <i>et al.</i> (1995) The kinetics and distribution of C9 and SC5b-9 <i>in vivo</i>: effects of complement activation. <a href="#">Clin Exp Immunol. 100 (1): 40-6.</a></li> <li>3. Rus, H. <i>et al.</i> (2005) C5b-9 complement complex in autoimmune demyelination and multiple sclerosis: dual role in neuroinflammation and neuroprotection. <a href="#">Ann Med. 37 (2): 97-104.</a></li> </ol>
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10334 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2611">https://www.bio-rad-antibodies.com/SDS/MCA2611</a> 10334
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

- Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
- Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
- Goat Anti Mouse IgG (STAR76...) [RPE](#)

Goat Anti Mouse IgG (STAR70...) [FITC](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),  
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),  
[FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
Goat Anti Mouse IgG (STAR77...) [HRP](#)

## Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M390493:210913'

Printed on 09 Dec 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)