

Datasheet: MCA2652

BATCH NUMBER 164778

Description:	MOUSE ANTI RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Specificity:	RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Other names:	RSV
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	0651
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/20 - 1/200
Immunoprecipitation			▪	
Western Blotting			▪	
Immunofluorescence	▪			1/10 - 1/50
Functional Assays			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Viral
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.1mg/ml
Immunogen	Fusion protein of the human Respiratory Syncytial Virus.
RRID	AB_844557
Specificity	<p>Mouse anti Respiratory syncytial virus antibody, clone 0651 detects the fusion protein of human respiratory syncytial virus (RSV), groups A and B. RSV is a negative-sense, single-stranded RNA virus and is a member of the <i>Paramyxoviridae</i> family. RSV causes respiratory tract infections in patients of all ages, but particularly affects infants and the immunosuppressed.</p> <p>RSV encodes three envelope glycoproteins, a small hydrophobic (SH) protein of unknown function, a glycoprotein (G) known as the attachment protein, and a fusion (F) protein. The F protein directs fusion of viral and cellular membranes, resulting in viral penetration, and can lead to the formation of syncytia.</p> <p>The F protein is thought to be the principal antigen responsible for inducing an immune response.</p>
References	<ol style="list-style-type: none"> Schlender, J. <i>et al.</i> (2003) Respiratory syncytial virus (RSV) fusion protein subunit F2, not attachment protein G, determines the specificity of RSV infection. J Virol. 77 (8): 4609-16. Kim, Y.K. <i>et al.</i> (2007) Genetic variability of the fusion protein and circulation patterns of genotypes of the respiratory syncytial virus. J Med Virol. 79 (6): 820-8.
Further Reading	<ol style="list-style-type: none"> Cianci, C. <i>et al.</i> (2005) Antiviral activity and molecular mechanism of an orally active respiratory syncytial virus fusion inhibitor. J Antimicrob Chemother. 55 (3): 289-92.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2652 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M403749:220803'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)