

## Datasheet: MCA2672GA

**BATCH NUMBER 161371**

<b>Description:</b>	MOUSE ANTI AMPK ALPHA 1/2
<b>Specificity:</b>	AMPK ALPHA 1/2
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	84.1
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Drosophila
<b>Species Cross Reactivity</b>	Based on sequence similarity, is expected to react with: Vertebrates, Human, Mouse, Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	29 amino acid peptide within the N-terminal catalytic domain of <i>Drosophila</i> AMPK alpha.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q13131</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">5562</a>    PRKAA1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	AMPK1
<b>RRID</b>	AB_1604624
<b>Specificity</b>	<p><b>Mouse anti AMPK Alpha 1/2 antibody, clone 84.1</b> recognises both the alpha 1 and the alpha 2 isoforms of AMP-activated protein kinase (AMPK). AMPK is a highly conserved heterotrimeric protein kinase found in all eukaryotes. The heterotrimeric complex exists as multiple isoforms and consists of a catalytic alpha subunit, which is stabilized by non-catalytic regulatory beta- and gamma- subunits. AMPK plays a vital role in maintaining cellular energy homeostasis and is activated by rising levels of AMP under conditions of hypoxia, ischemia, exercise, low glucose levels.</p> <p>AMPK may play an important role in the control of type 2 diabetes and metabolic syndromes. AMPK functions as a suppressor of cell proliferation of both malignant and normal cells (<a href="#">Motoshima et al. 2006</a>).</p>
<b>Western Blotting</b>	MCA2673GA detects a band of approximately 63kDa in rat brain cell lysates.
<b>Storage</b>	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2672GA">https://www.bio-rad-antibodies.com/SDS/MCA2672GA</a> 10040

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight@488</a> , <a href="#">DyLight@550</a> , <a href="#">DyLight@650</a> , <a href="#">DyLight@680</a> , <a href="#">DyLight@800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M383861:210513'

Printed on 18 Jan 2024