

Datasheet: MCA2672GA

Description:	MOUSE ANTI AMPK ALPHA 1/2
Specificity:	AMPK ALPHA 1/2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	84.1
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Drosophila

Species Cross Reactivity

Based on sequence similarity, is expected to react with: Vertebrates, Human, Mouse, Rat
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	29 amino acid peptide within the N-terminal catalytic domain of <i>Drosophila</i> AMPK alpha.
External Database Links	<p>UniProt: Q13131 Related reagents</p> <p>Entrez Gene: 5562 PRKAA1 Related reagents</p>
Synonyms	AMPK1
RRID	AB_1604624
Specificity	<p>Mouse anti AMPK Alpha 1/2 antibody, clone 84.1 recognises both the alpha 1 and the alpha 2 isoforms of AMP-activated protein kinase (AMPK). AMPK is a highly conserved heterotrimeric protein kinase found in all eukaryotes. The heterotrimeric complex exists as multiple isoforms and consists of a catalytic alpha subunit, which is stabilized by non-catalytic regulatory beta- and gamma- subunits. AMPK plays a vital role in maintaining cellular energy homeostasis and is activated by rising levels of AMP under conditions of hypoxia, ischemia, exercise, low glucose levels.</p> <p>AMPK may play an important role in the control of type 2 diabetes and metabolic syndromes. AMPK functions as a suppressor of cell proliferation of both malignant and normal cells (Motoshima et al. 2006).</p>
Western Blotting	MCA2673GA detects a band of approximately 63kDa in rat brain cell lysates.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2672GA 10040

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M405481:220916'

Printed on 25 Mar 2023
