

Datasheet: MCA2673GA

BATCH NUMBER 169037

Description:	MOUSE ANTI AMPK ALPHA 1/2
Specificity:	AMPK ALPHA 1/2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	34.2
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Drosophila

Species Cross Reactivity Reacts with: Human, Mouse, Rat
Based on sequence similarity, is expected to react with: Vertebrates
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form Purified IgG - liquid

Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	29 amino acid peptide within the N-terminal catalytic domain of Drosophila AMPK alpha.
External Database Links	<p>UniProt:</p> <p>P54646 Related reagents</p> <p>Q13131 Related reagents</p> <p>Entrez Gene:</p> <p>5563 PRKAA2 Related reagents</p> <p>5562 PRKAA1 Related reagents</p>
Synonyms	AMPK, AMPK1, AMPK2
RRID	AB_1604623
Specificity	<p>Mouse anti AMPK Alpha 1/2 antibody, clone 34.2 recognizes both the alpha 1 and the alpha 2 isoforms of AMP-activated protein kinase (AMPK). AMPK is a highly conserved heterotrimeric protein kinase found in all eukaryotes. The heterotrimeric complex exists as multiple isoforms and consists of a catalytic alpha subunit, which is stabilized by non-catalytic regulatory beta- and gamma- subunits. AMPK plays a vital role in maintaining cellular energy homeostasis and is activated by rising levels of AMP under conditions of hypoxia, ischemia, exercise, low glucose levels.</p> <p>AMPK may play an important role in the control of type 2 diabetes and metabolic syndromes. AMPK functions as a suppressor of cell proliferation of both malignant and normal cells (Motoshima et al. 2006).</p>
Western Blotting	MCA2673GA detects a band of approximately 63kDa in rat brain cell lysates.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2673GA 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M431588:240730'

Printed on 18 Mar 2025