

Datasheet: MCA2681PE

Description:	RAT ANTI MOUSE ABCA1:RPE
Specificity:	ABCA1
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	5A1-1422
Isotype:	IgG2a
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse						
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized						
Reconstitution	Reconstitute with 1ml distilled water						
Max Ex/Em	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>RPE 488nm laser</td> <td>496</td> <td>578</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	RPE 488nm laser	496	578
Fluorophore	Excitation Max (nm)	Emission Max (nm)					
RPE 488nm laser	496	578					
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant						
Buffer Solution	Phosphate buffered saline						
Preservative	0.09% Sodium Azide (NaN ₃)						
Stabilisers	1% Bovine Serum Albumin 5% Sucrose						

Immunogen	ABCA1 transfected HeLa cells.
External Database Links	<p>UniProt: P41233 Related reagents</p> <p>Entrez Gene: 11303 Abca1 Related reagents</p>
Synonyms	Abc1
RRID	AB_2220006
Specificity	<p>Rat anti Mouse ABCA1 antibody, clone 5A1.1422 recognizes murine adenosine triphosphate (ATP) Binding cassette transporter 1 (ABCA1). The ABC transporters are a large family of conserved proteins that transport a wide variety of molecules across cellular membranes. ABCA1 is a member of the ABC-A sub-family, which acts as a lipid translocator. The molecule was originally identified as a scavenger receptor on macrophages and research shows that ABCA1 also plays a major role in cholesterol metabolism. ABCA1 may play an important role in protecting against cardiovascular disease.</p> <p>Mutations in ABCA1 gene have been associated with Tangier disease, a genetic disorder of lipid metabolism, and the more common Familial HDL deficiency, presenting with many of the same symptoms.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Pradel, L.C. <i>et al.</i> (2009) ATP-binding cassette transporter hallmarks tissue macrophages and modulates cytokine-triggered polarization programs. Eur J Immunol. 39 (8): 2270-80. 2. Wang, Y. <i>et al.</i> (2019) Smooth Muscle Cells Contribute the Majority of Foam Cells in ApoE (Apolipoprotein E)-Deficient Mouse Atherosclerosis. Arterioscler Thromb Vasc Biol. : ATVBAHA119312434. [Epub ahead of print]
Further Reading	1. Zarubica, A. <i>et al.</i> (2007) ABCA1, from pathology to membrane function. Pflugers Arch. 453 (5): 569-79.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20487 available at: 20487: https://www.bio-rad-antibodies.com/uploads/MSDS/20487.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:RPE \(MCA1212PE\)](#)

Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)

[MOUSE SEROBLOCK FcR \(BUF041B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M375524:210104'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)