

Datasheet: MCA2694

| | |
|----------------------|--------------------------|
| Description: | RAT ANTI MOUSE CD8 ALPHA |
| Specificity: | CD8 ALPHA |
| Other names: | LY-2 |
| Format: | Purified |
| Product Type: | Monoclonal Antibody |
| Clone: | 53-6.7 |
| Isotype: | IgG2a |
| Quantity: | 0.25 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | 1/20 |
| Immunohistology - Frozen | ▪ | | | |
| Immunohistology - Paraffin | | | ▪ | |
| ELISA | | | ▪ | |
| Immunoprecipitation | ▪ | | | |
| Western Blotting | | | ▪ | |
| Immunofluorescence | ▪ | | | |
| Functional Assays | | | ▪ | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|---|
| Target Species | Mouse |
| Product Form | Purified IgG - liquid |
| Preparation | Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | <0.1% Sodium Azide (NaN ₃) |

| | |
|---------------------------------------|--|
| Approx. Protein Concentrations | IgG concentration 0.5mg/ml |
| Immunogen | Mouse thymus or spleen. |
| External Database Links | <p>UniProt: P01731 Related reagents</p> <p>Entrez Gene: 12525 Cd8a Related reagents</p> |
| Synonyms | Lyt2, Lyt-2 |
| RRID | AB_905993 |
| Specificity | <p>Rat anti Mouse CD8 alpha antibody, clone 53-6.7 detects mouse CD8a, also known as Lyt-2. CD8α is the alpha chain of the CD8 antigen. CD8 is expressed on T-cells, and exists either a homodimer composed of two alpha chains, or a heterodimer composed of one alpha and one beta chain.</p> <p>CD8 functions as a co-receptor with the T-cell receptor and mediates efficient cell to cell interactions within the immune system. The CD8 alpha chain specifically binds to class I MHC molecules.</p> <p>Rat anti Mouse CD8 alpha antibody, clone 53-6.7 has been reported to block antigen presentation via MHC class I and inhibit T cell responses to IL-2. Rat anti Mouse CD8 alpha antibody, clone 53-6.7 has also been used for depletion of CD8a⁺ cells.</p> |
| Flow Cytometry | <p>Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.</p> <p>The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors. This may be reduced by using SeroBlock FcR (BUF041A/B).</p> |
| Immunohistology | MCA2694 has been reported as working well on Zinc-fixed paraffin-embedded sections. However, inconsistent results have been reported on formalin-fixed paraffin-embedded tissue sections. |
| References | <ol style="list-style-type: none"> 1. Takahashi, K. <i>et al.</i> (1992) CD4 and CD8 regulate interleukin 2 responses of T cells. Proc Natl Acad Sci U S A. 89 (12): 5557-61. 2. Shore, D.A. <i>et al.</i> (2008) The crystal structure of CD8 in complex with YTS156.7.7 Fab and interaction with other CD8 antibodies define the binding mode of CD8 alphabeta to MHC class I. J Mol Biol. 384 (5): 1190-202. 3. McDole, J.R. <i>et al.</i> (2010) Rapid formation of extended processes and engagement of Theiler's virus-infected neurons by CNS-infiltrating CD8 T cells. Am J Pathol. 177 (4): 1823-33. 4. Schuster, C. <i>et al.</i> (2013) A missing PD-L1/PD-1 coinhibition regulates diabetes induction by preproinsulin-specific CD8 T-cells in an epitope-specific manner. PLoS One. |

[8 \(8\): e71746.](#)

5. Johnson, H.L. *et al.* (2014) Perforin competent CD8 T cells are sufficient to cause immune-mediated blood-brain barrier disruption. [PLoS One. 9 \(10\): e111401.](#)

6. Zhang, M.Z. *et al.* (2015) Inhibition of cyclooxygenase-2 in hematopoietic cells results in salt-sensitive hypertension. [J Clin Invest. 125 \(11\): 4281-94.](#)

Further Reading 1. Jabbari, A. & Harty, J.T. (2006) The generation and modulation of antigen-specific memory CD8 T cell responses. [J Leukoc Biol. 80 \(1\): 16-23.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2694>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

| | |
|--|---|
| Goat Anti Rat IgG (STAR73...) | RPE |
| Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...) | DyLight@550 , DyLight@650 , DyLight@800 |
| Rabbit Anti Rat IgG (STAR21...) | HRP |
| Rabbit Anti Rat IgG (STAR16...) | DyLight@800 |
| Goat Anti Rat IgG (STAR131...) | Alk. Phos. , Biotin |
| Rabbit Anti Rat IgG (STAR17...) | FITC |
| Goat Anti Rat IgG (STAR72...) | HRP |
| Goat Anti Rat IgG (STAR69...) | FITC |

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL \(MCA1212\)](#)

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M390159:210825'

Printed on 13 Sep 2023