

Datasheet: MCA2712B

Description:	MOUSE ANTI HUMAN ANNEXIN V:Biotin
Specificity:	ANNEXIN V
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	VAA-33
Isotype:	IgG2a
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Immunogen	Recombinant Annexin V.

External Database Links

UniProt:

[P08758](#)

[Related reagents](#)

Entrez Gene:

Synonyms	ANX5, ENX2, PP4
-----------------	-----------------

RRID	AB_915142
-------------	-----------

Specificity	<p>Mouse anti Human Annexin V antibody, clone VAA-33 recognizes human Annexin V, a phospholipid binding protein which binds with high affinity and selectivity to phosphatidylserine (PS), in the presence of calcium.</p> <p>The translocation of PS from the cytosolic leaflet to the outer leaflet of the plasma membrane, due to a loss in plasma membrane asymmetry during apoptosis, acts as an in vivo signal for the removal of the dying cell by phagocytosis. Detection of PS using Annexin V therefore provides a reliable diagnostic tool for the indication of early programmed cell death.</p> <p>Annexin V has also been identified as an anticoagulant protein in the blood coagulation cascade, by acting as an inhibitor of prothrombin activation. The presence of antibodies to Annexin V, are associated with systemic lupus erythematosus (SLE), recurrent spontaneous abortions and systemic sclerosis (SSc).</p> <p>Mouse anti Human Annexin V antibody, clone VAA-33 has been reported to work in Western blotting applications (Wedher et al. 2009).</p>
--------------------	---

Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.
-----------------------	---

Further Reading	1. Gerke, V. & Moss, S.E. (2002) Annexins: from structure to function. Physiol Rev. 82 (2): 331-71.
------------------------	---

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10058 available at: 10058: https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[ANNEXIN V:FITC ASSAY KIT \(ANNEX100F\)](#)

[MOUSE ANTI HUMAN ANNEXIN V \(MCA2712\)](#)

[MOUSE ANTI HUMAN ANNEXIN V:FITC \(MCA2712F\)](#)

[ANNEXIN V:FITC ASSAY KIT \(ANNEX300F\)](#)

[ANNEXIN V:Biotin ASSAY KIT \(ANNEX100B\)](#)

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.

Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M367359:200529'

Printed on 01 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)