

Datasheet: MCA2747GA

BATCH NUMBER 153274

Description:	RAT ANTI MOUSE CD83
Specificity:	CD83
Other names:	HB15
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	3D11
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)

Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
External Database Links	<p>UniProt: O88324 Related reagents</p> <p>Entrez Gene: 12522 Cd83 Related reagents</p>
RRID	AB_10614522
Fusion Partners	Spleen cells from immunised rats were fused with cells of the P3X63.Ag8.653 mouse myeloma cell line.
Specificity	Rat anti Mouse CD83 antibody, clone 3D11 recognizes murine CD83, a 45 kDa type 1 glycoprotein belonging to the immunoglobulin superfamily. Originally described as a maturation marker for dendritic cells, CD83 plays a role in the regulation of antigen presentation, thymic T-cell maturation, B-cell function and peripheral T-cell activation. Activated dendritic cells and B-cells release a soluble form of CD83 that is reported to have an immunosuppressive effect, suggesting a therapeutic role for sCD83 in the regulation of immune disorders.
Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.
References	<ol style="list-style-type: none"> 1. Fujimoto, Y. & Tedder, T.F. (2006) CD83: a regulatory molecule of the immune system with great potential for therapeutic application. J Med Dent Sci. 53 (2): 85-91. 2. He, Y.T. <i>et al.</i> (2016) <i>In vitro</i> generation of cytotoxic T lymphocyte response using dendritic cell immunotherapy in osteosarcoma. Oncol Lett. 12 (2): 1101-6. 3. El-Ashmawy, N.E. <i>et al.</i> (2015) <i>In vitro</i> and <i>in vivo</i> studies of the immunomodulatory effect of <i>Echinacea purpurea</i> on dendritic cells. Journal of Genetic Engineering and Biotechnology. 13 (2): 185-92.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2747GA 10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...)	DyLight®800
Rabbit Anti Rat IgG (STAR17...)	FITC
Goat Anti Rat IgG (STAR72...)	HRP
Goat Anti Rat IgG (STAR69...)	FITC
Goat Anti Rat IgG (STAR73...)	RPE
Rabbit Anti Rat IgG (STAR21...)	HRP
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	DyLight®550 , DyLight®650 , DyLight®800
Goat Anti Rat IgG (STAR131...)	Alk. Phos. , Biotin

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL \(MCA1212\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M367395:200529'

Printed on 18 Jan 2024