

Datasheet: MCA278

BATCH NUMBER 155847

Description:	MOUSE ANTI RAT IgG2a HEAVY CHAIN
Specificity:	IgG2a HEAVY CHAIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	MARG2a-1
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1.25ug/ml as coating antibody
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide
Approx. Protein	IgG concentration 1 mg/ml

Concentrations

Immunogen purified rat IgG2a.

External Database Links

UniProt:

[P20760](#) [Related reagents](#)

Entrez Gene:

[679045](#) LOC679045 [Related reagents](#)

RRID AB_322914

Fusion Partners Spleen cells from immunised BALB/c mice were fused with cells of the PAI-0 mouse myeloma cell line.

Specificity **Mouse anti Rat IgG2a heavy chain antibody, clone MARG2a-1** recognizes the rat gamma 2a immunoglobulin heavy chain, and does not cross-react with other immunoglobulin classes or subclasses.

Avidity of Mouse anti Rat IgG2a heavy chain antibody, clone MARG2a-1 for rat IgG2a = $2.1 \times 10^9 \text{ M}^{-1}$.

ELISA

This product may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with [MCA1296P](#) as detection reagent.

Purified rat IgA ([PRP04](#)) may be used as a standard.

References

1. Bazin, H. *et al.* (1974) Three classes and four (sub)classes of rat immunoglobulins: IgM, IgA, IgE and IgG1, IgG2a, IgG2b, IgG2c. [Eur J Immunol. 4 \(1\): 44-8.](#)
 2. Shulman, M. *et al.* (1978) A better cell line for making hybridomas secreting specific antibodies. [Nature. 276 \(5685\): 269-70.](#)
 3. Bazin, H. *et al.* (1973) Brief communication: Transplantable immunoglobulin-secreting tumors in rats. V. Monoclonal immunoglobulins secreted by 250 ileocecal immunocytomas in LOU-Wsl rats. [J Natl Cancer Inst. 51 \(4\): 1359-61.](#)
 4. Sandin, L.C. *et al.* (2014) Locally delivered CD40 agonist antibody accumulates in secondary lymphoid organs and eradicates experimental disseminated bladder cancer. [Cancer Immunol Res. 2 \(1\): 80-90.](#)
 5. Bézie S *et al.* (2015) Compensatory Regulatory Networks between CD8 T, B, and Myeloid Cells in Organ Transplantation Tolerance. [J Immunol. 195 \(12\): 5805-15.](#)
 6. Bézie, S. *et al.* (2015) Fibrinogen-like protein 2/fibroleukin induces long-term allograft survival in a rat model through regulatory B cells. [PLoS One. 10 \(3\): e0119686.](#)
-

Further Reading

1. Köhler, G. & Milstein, C. (1975) Continuous cultures of fused cells secreting antibody of predefined specificity. [Nature. 256 \(5517\): 495-7.](#)
 2. Bazin, H. *et al.* (1978) Transplantable IgD immunoglobulin-secreting tumors in rat. [J. Immunol. 121: 2077-2082.](#)
-

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA278>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M367451:200529'

Printed on 18 Apr 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)