

Datasheet: MCA2809PE

Description:	MOUSE ANTI HUMAN CD5:RPE
Specificity:	CD5
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	UCHT2
Isotype:	IgG1
Quantity:	100 TESTS/0.5ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Chimpanzee, Cynomolgus monkey, Rhesus Monkey, Marmoset
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578

Preparation

Purified IgG prepared by affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

External Database

Links

UniProt:

[P06127](#) [Related reagents](#)

Entrez Gene:

[921](#) CD5 [Related reagents](#)

Synonyms

LEU1

Specificity

Mouse anti Human CD5 antibody, clone UCHT2 recognizes the human CD5 cell surface antigen. CD5 is expressed by T cells and by a small subset of B cells. CD5 is also expressed by certain types of B cell malignancy, including B cell chronic lymphocytic leukaemia (B-CLL).

Flow Cytometry

Use 5ul of the suggested working dilution to label 10⁶ cells or 100ul whole blood

References

1. Renaudineau, Y. *et al.* (2005) An alternative exon 1 of the CD5 gene regulates CD5 expression in human B lymphocytes. [Blood 106 \(8\): 2781-9.](#)

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M374945:201209'

Printed on 22 Mar 2021
