

Datasheet: MCA2818

**BATCH NUMBER 180626**

<b>Description:</b>	MOUSE ANTI HUMAN PAPILOMAVIRUS 16 ONCOPROTEIN E7
<b>Specificity:</b>	PAPILOMAVIRUS 16 ONCOPROTEIN E7
<b>Other names:</b>	HPV
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	716-D1
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Viral
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium azide

<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	Oncoprotein E7 from human papillomavirus 16.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P03129</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">1489079</a> E7    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_1102614
<b>Specificity</b>	<p><b>Mouse anti Human Papillomavirus 16 Oncoprotein E7 antibody, clone 716-D1</b> recognizes both the monomer and dimer forms of the E7 oncoprotein of human papilloma virus 16. Human papillomavirus (HPV) is a diverse group of DNA-based viruses that infect skin and mucous membranes of humans and animals. Some HPV types are the causative agents of cervical cancer, with types 16 and 18 being particularly high-risk. The viral proteins E6 and E7 disrupts normal cell cycle regulation by interacting with p53 (a tumor-suppressing transcription factor) and Rb (retinoblastoma protein, also a tumor-suppressor). E7 particularly binds to Rb and histone deacetylases, resulting in activation of the E2F genes, which code for a family of transcription factors. The viral proteins E6 and E7 might be of particular interest in the development of therapeutic vaccines, since they are expressed early in viral infection.</p> <p>Mouse anti Human Papillomavirus 16 Oncoprotein E7 antibody, clone 716-D1 cross-reacts with both the monomer and dimer forms of the E7 protein of HPV-18.</p>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10040 available at:  <a href="https://www.bio-rad-antibodies.com/SDS/MCA2818">https://www.bio-rad-antibodies.com/SDS/MCA2818</a>  10040</p>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>

Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)  
Goat Anti Mouse IgG (STAR76...) [RPE](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),  
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),  
[FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M367484:200529'

**Printed on 25 Mar 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)