

Datasheet: MCA2828

**BATCH NUMBER 155439**

<b>Description:</b>	MOUSE ANTI RABIES VIRUS GLYCOPROTEIN
<b>Specificity:</b>	RABIES VIRUS GLYCOPROTEIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	1C5
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.2 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	
Immunofluorescence	▪			
Functional Assays (1)	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.  
(1)Sodium azide needs to be removed by dialysis before use of this antibody in functional assay

<b>Target Species</b>	Viral
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> )

<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml
<b>Immunogen</b>	Purified rabies virus.
<b>RRID</b>	AB_1125351
<b>Fusion Partners</b>	Spleen cells from immunised Balb/c mice were fused with cells of the Sp2/0 myeloma cells.
<b>Specificity</b>	<p><b>Mouse anti Rabies Virus Glycoprotein antibody, clone 1C5</b> recognizes the glycoprotein of rabies virus, a Lyssavirus that causes rabies in mammals. The virus is found in the saliva of infected animals and is mostly spread through bites. The virus then travels through the peripheral nervous system to the central nervous system. Once it reaches the brain it causes acute encephalitis. Symptoms include paralysis, cerebral dysfunction, anxiety, insomnia, confusion, anger, paranoia, terror and hallucination. The production of a large amount of saliva, the inability to speak or swallow and associated hydrophobia is also common in the later stages of the disease. The time between exposure to the virus and onset of symptoms can range from about 2 weeks to several months. Once symptoms appear, death is virtually inevitable. However, the disease can be effectively treated by post-exposure prophylaxis, through vaccination, provided this is done as soon as possible.</p> <p>Mouse anti Rabies Virus Glycoprotein antibody, clone 1C5 recognizes the purified glycoprotein of the rabies virus as well as the native virus. It has been shown to react with intact infected BHK cells and has therapeutic effect on infected mice. The antibody has virus-neutralising activity with several rabies virus strains such as CVS, Lagos bat, Mokola, Flury LEP and Duvenhage.</p>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2828">https://www.bio-rad-antibodies.com/SDS/MCA2828</a></p> <p>10040</p>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)  
 Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)  
 Goat Anti Mouse IgG (STAR76...) [RPE](#)  
 Goat Anti Mouse IgG (STAR70...) [FITC](#)  
 Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)  
 Goat Anti Mouse IgG (STAR77...) [HRP](#)  
 Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
 Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
 Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
 'M378860:210309'

**Printed on 26 Jun 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)