

Datasheet: MCA2831

BATCH NUMBER 150425

Description:	MOUSE ANTI SALMONELLA PARATYPHI A LPS GROUP A 0-2 ANTIGEN
Specificity:	SALMONELLA PARATYPHI A LPS GROUP A 0-2 ANTIGEN
Other names:	O-ANTIGEN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	10B10G
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting			▪	
Immunofluorescence	▪			
Immunoblotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide (NaN ₃)

Stabilisers

Approx. Protein Concentrations IgG concentration 1.0mg/ml

Immunogen Heat inactivated cells of *Salmonella paratyphi A*.

RRID AB_1125359

Fusion Partners Spleen cells from immunised Balb/c mice were fused with cells of Sp2/0 the myeloma cells.

Specificity **Mouse anti Salmonella paratyphi A LPS Group A 0-2 Antigen antibody, clone 10B10G** recognizes the O antigen of *Salmonella enterica*, serovar paratyphi A (also known as *S. paratyphi*). The O antigen is a polysaccharide that extends from the core oligosaccharide, which itself is attached to lipid A. These three components together form the lipopolysaccharide (LPS), a major component of the outer membrane of Gram-negative bacteria. The LPS is an endotoxin and induces a strong immune response from the host. Its mutation or removal results in the death of the bacterial cell and as such is a prime target for antimicrobial drugs. The O antigen is the outermost component of the LPS and functions as a receptor for bacteriophages. It is also important in the immune response of the host. Considerable diversity is seen within *Salmonella* O antigens.

Salmonella enterica, serovar paratyphi A belongs to serogroup A of Salmonella and causes paratyphoid fever, which is similar to typhoid fever but more benign.

Mouse anti Salmonella paratyphi A LPS Group A 0-2 Antigen antibody, clone 10B10G recognizes the 0-2 antigenic determinant of *Salmonella* serogroup A and does not cross-react with the LPS of serogroups B, D or E of *Salmonella* spp.

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2831>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...) [HRP](#)

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M367491:200529'

Printed on 27 Apr 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)