

Datasheet: MCA3263Z

BATCH NUMBER 144519

Description:	MOUSE ANTI HUMAN MIF:Preservative Free
Specificity:	MIF
Format:	Preservative Free
Product Type:	Monoclonal Antibody
Clone:	2A10-4D3
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			0.1 - 10 ug/ml
ELISA			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant protein with a GST tag corresponding to amino acids 1 - 116 of human MIF.
External Database Links	<p>UniProt: P14174 Related reagents</p> <p>Entrez Gene: 4282 MIF Related reagents</p>
Synonyms	GLIF, MMIF
RRID	AB_2226469
Fusion Partners	Spleen cells from BALB/c mice were fused with cells from the Sp2/0 myeloma cell line.
Specificity	<p>Mouse anti Human MIF antibody, clone 2A10-4D3 recognizes human Macrophage migration inhibitory factor, also known as MIF, glycosylation-inhibiting factor, L-dopachrome isomerase or phenylpyruvate tautomerase. MIF is a 115 amino acid ~12.5 kDa pro-inflammatory cytokine playing a role in the innate immune response to bacterial pathogens. MIF plays a role in the regulation of macrophage function in host defense through the suppression of anti-inflammatory effects of glucocorticoids (UniProt:14174). MIF and the JAB1 protein form a complex in the cytosol near the peripheral plasma membrane, which may indicate an additional role in integrin signaling pathways (Kleemann et al. 2000). Mutation of the MIF gene is implicated in susceptibility and severity of early onset rheumatoid arthritis (Llamas-Covarrubias <i>et al.</i> 2013).</p> <p>Mouse anti Human MIF antibody, clone 2A10-4D3 has been utilized successfully for the micro immunoaffinity isolation of human MIF from human serum for subsequent mass spectrometric immunoassay (Sherma 2014).</p>
References	1. Sherma, N.D. <i>et al.</i> (2014) Mass Spectrometric Immunoassay for the qualitative and quantitative analysis of the cytokine Macrophage Migration Inhibitory Factor (MIF). Proteome Sci. 12 (1): 52.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/MCA3263Z 10162

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M367572:200529'

Printed on 25 Mar 2023
