

## Datasheet: MCA403

**BATCH NUMBER 159658**

<b>Description:</b>	MOUSE ANTI INFLUENZA B NUCLEOPROTEIN
<b>Specificity:</b>	INFLUENZA B NUCLEOPROTEIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	B017 (B35G)
<b>Isotype:</b>	IgG2b
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Viral
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography from tissue culture supernatant.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein</b>	IgG concentration 1.0 mg/ml

## Concentrations

**Immunogen** Influenza B/Lee/40 and B/Singapore/-222/79 viruses.

**RRID** AB\_2298473

**Fusion Partners** Spleen cells from immunised BALB/c mice were fused with cells of the P3 Ag8.653 mouse myeloma cell line.

**Specificity** **Mouse anti Influenza B Nucleoprotein antibody, clone B017** recognises an epitope within the nucleoprotein of influenza B virus.

The influenza viruses, classified as type A, B and C, are members of the *Orthomyxoviridae* family which differ in their epidemiology and host ranges and lack serological cross-reactivity of their internal components, especially their matrix proteins and nucleoprotein. Influenza B virus is a slow-mutating single stranded RNA virus subject to antigenic drift which, although enough to prevent lasting immunity, prevents influenza B from causing pandemics.

Mouse anti Influenza B Nucleoprotein antibody, clone B017 can be used in influenza B IFA typing in conjunction with [MCA2717](#) (clone B114).

## References

1. Walls, H.H. *et al.* (1986) Characterization and evaluation of monoclonal antibodies developed for typing influenza A and influenza B viruses. [J Clin Microbiol. 23 \(2\): 240-5.](#)
2. Zhirnov, O.P. *et al.* (1999) Caspase-dependent N-terminal cleavage of influenza virus nucleocapsid protein in infected cells. [J Virol. 73 \(12\): 10158-63.](#)
3. Ehrhardt, C. *et al.* (2007) Activation of phosphatidylinositol 3-kinase signaling by the nonstructural NS1 protein is not conserved among type A and B influenza viruses. [J Virol. 81 \(21\): 12097-100.](#)
4. Goujon, C. & Malim, M.H. (2010) Characterization of the alpha interferon-induced postentry block to HIV-1 infection in primary human macrophages and T cells. [J Virol. 84 \(18\): 9254-66.](#)
5. Dauber, B. *et al.* (2009) Influenza B virus ribonucleoprotein is a potent activator of the antiviral kinase PKR. [PLoS Pathog. 5 \(6\): e1000473.](#)
6. Lesch, M. *et al.* (2019) RNAi-based small molecule repositioning reveals clinically approved urea-based kinase inhibitors as broadly active antivirals [PLoS Pathog. 15\(3\):e1007601.](#)

## Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Guarantee** 12 months from date of despatch

**Health And Safety** Material Safety Datasheet documentation #10040 available at:

**Information** <https://www.bio-rad-antibodies.com/SDS/MCA403>  
10040

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)  
Goat Anti Mouse IgG (STAR76...) [RPE](#)  
Goat Anti Mouse IgG (STAR70...) [FITC](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
Goat Anti Mouse IgG (STAR77...) [HRP](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M389522:210806'

Printed on 24 Feb 2025

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)