

Datasheet: MCA423

BATCH NUMBER 180613

Description:	RAT ANTI MOUSE IgG3 HEAVY CHAIN
Specificity:	IgG3 HEAVY CHAIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	LO-MG3-7
Isotype:	IgM
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			5ug/ml as a coating antibody
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgM - liquid
Preparation	Purified IgM prepared by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide
Approx. Protein	IgM concentration 1 mg/ml

Concentrations

Immunogen Polyclonal and monoclonal mouse IgG3.

External Database Links

UniProt:

[P03987](#) [Related reagents](#)

Entrez Gene:

[380795](#) A1324046 [Related reagents](#)

RRID AB_321849

Specificity **Rat anti Mouse IgG3 Heavy Chain antibody, clone LO-MG3-7** recognizes the gamma 3 heavy chain of mouse IgG and does not cross react with other subclasses of murine IgG.

References

1. Bazin, H. (1982) Production of rat monoclonal antibodies with the LOU rat non secreting IR983F myeloma cell line. Prot Biol Fluids. Peters Ed: 29th Colloquim (1981). Pergamon Press, Oxford and N.Y., pp 615-8.
 2. Bazin, H. *et al.* (1984) Rat monoclonal antibodies. I. Rapid purification from *in vitro* culture supernatants. [J Immunol Methods. 66 \(2\): 261-9.](#)
 3. Bazin, H. *et al.* (1984) Rat monoclonal antibodies. II. A rapid and efficient method of purification from ascitic fluid or serum. [J Immunol Methods. 71 \(1\): 9-16.](#)
 4. Zivkovic, I. *et al.* (2012) Induction of decreased fecundity by tetanus toxoid hyper-immunization in C57BL/6 mice depends on the applied adjuvant. [Innate Immun. 18: 333-42.](#)
 5. DeGiovanni, C. *et al.* (2019) Cancer Vaccines Co-Targeting HER2/Neu and IGF1R. [Cancers \(Basel\). 11 \(4\)Apr 11 \[Epub ahead of print\].](#)
-

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA423>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rat IgM (302001...) [Biotin](#)

Goat Anti Rat IgM (STAR116...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M367681:200529'

Printed on 20 Oct 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)