

Datasheet: MCA4265Z

**BATCH NUMBER 164708**

<b>Description:</b>	MOUSE ANTI HUMAN SQSTM1:Preservative Free
<b>Specificity:</b>	SQSTM1
<b>Other names:</b>	p62
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	2C11
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			
Immunofluorescence	▪			0.1 - 10 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from ascites
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Ig concentration 0.5 mg/ml
<b>Immunogen</b>	Recombinant protein corresponding to aa 1 - 441 of human SQSTM1.

**External Database  
Links**

**UniProt:**

[Q13501](#)   [Related reagents](#)

**Entrez Gene:**

[8878](#)   SQSTM1   [Related reagents](#)

---

**Synonyms**

ORCA, OSIL

---

**RRID**

AB\_2286694

---

**Fusion Partners**

Spleen cells from BALB/c mice were fused with cells from the Sp2/0 myeloma cell line.

---

**Specificity**

**Mouse anti Human SQSTM1 antibody, clone 2C11** recognizes human SQSTM1, also known as Sequestosome-1, p62 or EBIAP. SQSTM1 is a ~62kDa ubiquitously expressed molecule bearing an N-terminal [OPR](#) domain which mediates homo-oligomerization and interactions with other proteins. SQSTM1 also bears a C-terminal [UBA](#) domain, responsible for binding to poly-ubiquitinated substrates.

SQSTM1 is involved in the autophagic degradation of poly-ubiquitinated proteins ([Wei et al. 2013](#)) and in bone remodelling through activation of osteoclasts. Mutations in the SQSTM1 gene have been implicated in Paget's disease of bone ([Laurin et al. 2002](#))

Mouse anti human SQSTM1, clone 2C11 has been used successfully for the identification of SQSTM1 in Hela cells by western blotting.

---

**Storage**

Store at -20°C only.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10162 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA4265Z>  
10162

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)   [HRP](#)  
Rabbit Anti Mouse IgG (STAR12...)   [RPE](#)  
Goat Anti Mouse IgG (STAR70...)   [FITC](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...)   [FITC](#)

Goat Anti Mouse IgG (STAR76...)

[RPE](#)

Goat Anti Mouse IgG (H/L) (STAR117...)

[Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR13...)

[HRP](#)

Goat Anti Mouse IgG (Fc) (STAR120...)

[FITC](#), [HRP](#)

## Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

**North & South** Tel: +1 800 265 7376

**Worldwide**

Tel: +44 (0)1865 852 700

**Europe**

Tel: +49 (0) 89 8090 95 21

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M405936:220916'

Printed on 25 Mar 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)