

Datasheet: MCA4340Z

BATCH NUMBER 173783

Description:	MOUSE ANTI HUMAN JARID1B:Preservative Free
Specificity:	JARID1B
Other names:	LYSINE SPECIFIC DEMETHYLASE 5B
Format:	Preservative Free
Product Type:	Monoclonal Antibody
Clone:	1G10
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin (1)	▪			0.1 - 10 ug/ml
Western Blotting	▪			0.1 - 10 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	<p>Reacts with: Rat, Mouse</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from ascites

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Ig concentration 0.5 mg/ml
Immunogen	Recombinant protein corresponding to aa 231-320 of human JARID1B
External Database Links	<p>UniProt: Q9UGL1 Related reagents</p> <p>Entrez Gene: 10765 KDM5B Related reagents</p>
Synonyms	JARID1B, PLU1, RBBP2H1
RRID	AB_1221667
Fusion Partners	Spleen cells from BALB/c mice were fused with cells from the Sp2/0 myeloma cell line.
Specificity	Mouse anti Human JARID1B antibody, clone 1G10 recognizes human Lysine-specific demethylase 5B, also known as JARID1B, Cancer/testis antigen 31, Histone demethylase JARID1B, PLU-1, Retinoblastoma-binding protein 2 homolog 1 or Jumonji/ARID domain-containing protein 1B. JARID1B is a 1544 amino acid ~180 kDa ubiquitously expressed histone demethylase.
References	<ol style="list-style-type: none"> 1. Kuźbicki, L. <i>et al.</i> (2013) JARID1B expression in human melanoma and benign melanocytic skin lesions. Melanoma Res. 23 (1): 8-12. 2. Kuźbicki, Ł. <i>et al.</i> (2017) Prognostic significance of RBP2-H1 variant of JARID1B in melanoma. BMC Cancer. 17 (1): 854. 3. Di Nisio, E. <i>et al.</i> (2023) A truncated and catalytically inactive isoform of KDM5B histone demethylase accumulates in breast cancer cells and regulates H3K4 tri-methylation and gene expression. Cancer Gene Ther. Jan 26 [Epub ahead of print].
Further Reading	1. Xiang, Y. <i>et al.</i> (2007) JARID1B is a histone H3 lysine 4 demethylase up-regulated in prostate cancer. Proc Natl Acad Sci U S A. 104: 19226-31.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety	Material Safety Datasheet documentation #10162 available at:

Information <https://www.bio-rad-antibodies.com/SDS/MCA4340Z>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M405939:220916'

Printed on 26 May 2026