

Datasheet: MCA4646SBUV700

Description:	MOUSE ANTI HUMAN CD335:StarBright UltraViolet 700
Specificity:	CD335
Other names:	NKp46
Format:	StarBright UltraViolet 700
Product Type:	Monoclonal Antibody
Clone:	9E2
Isotype:	IgG1
Quantity:	100 TESTS/0.5ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Staining at lower temperatures generally requires longer incubation to reach antibody binding equilibrium, and this effect is more pronounced for StarBright Dyes due to their larger size. For optimum staining, Bio-Rad recommends incubating samples either overnight at 4°C or 1 hour at room temperature.

Target Species	Human		
Product Form	Purified IgG conjugated to StarBright UltraViolet 700 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	StarBright UltraViolet 700	340	700
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		

Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 0.1% Pluronic F68 0.1% PEG 3350 0.05% Tween 20
Approx. Protein Concentrations	For information on the concentration of our StarBright Dye conjugated reagents please visit our FAQ page.
Immunogen	CD335-Fc fusion protein
External Database Links	UniProt: O76036 Related reagents Entrez Gene: 9437 NCR1 Related reagents
Synonyms	LY94
Specificity	Mouse anti Human CD335 antibody, clone 9E2 recognizes human NK (natural killer) cell-specific marker CD335, otherwise known as NKp46, a 46kDa type I transmembrane glycoprotein, and the first identified member of the natural cytotoxicity receptor (NCR) family. CD335 is uniquely expressed by all resting and activated NK cells, and considered to be the major NK cell lysis receptor for autologous pathogen-infected and tumour target cells, during natural cytotoxicity responses. Ligands of CD335 include viral haemagglutinins (HAs), and heparan sulphate proteoglycans (HSPGs) on the surface of tumour cells. Mouse anti Human CD335 antibody, clone 9E2 mediates re-directed NK cell cytotoxicity (Yusa <i>et al.</i> 2002, Chen <i>et al.</i> 2007).
Flow Cytometry	Use 5µl of the suggested working dilution to label 0.5x10 ⁶ cells in 100µl. Best practices suggest a 5 min centrifugation at 6,000g prior to sample application.
References	1. Yusa, S. <i>et al.</i> (2002) SHP-1- and phosphotyrosine-independent inhibitory signaling by a killer cell Ig-like receptor cytoplasmic domain in human NK cells. J Immunol. 168 (10): 5047-57. 2. Chklovskaya, E. <i>et al.</i> (2004) Reconstitution of dendritic and natural killer-cell subsets after allogeneic stem cell transplantation: effects of endogenous flt3 ligand. Blood. 103 (10): 3860-8. 3. Chen, Y. <i>et al.</i> (2007) Prostaglandin D2 suppresses human NK cell function via signaling through D prostanoid receptor. J Immunol. 179 (5): 2766-73.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.

Guarantee	12 months from date of despatch
Acknowledgements	This product is covered by U.S. Patent No. 10,150,841 and related U.S. and foreign counterparts
Health And Safety Information	Material Safety Datasheet documentation #20471 available at: https://www.bio-rad-antibodies.com/SDS/MCA4646SBUV700
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M452395:260710'

Printed on 10 Jul 2026