

Datasheet: MCA4684

BATCH NUMBER 0709

Description:	MOUSE ANTI BOVINE GFAP
Specificity:	GFAP
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	4A11
Isotype:	IgG2b
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/500 - 1/10000
Immunofluorescence			▪	
Radioimmunoassays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Human, Mouse, Rat, Sheep, Dog, Pig, Rabbit, Guinea Pig, Chicken
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Bovine spinal cord homogenate.
External Database Links	<p>UniProt:</p> <p>P14136 Related reagents</p> <p>P47819 Related reagents</p> <p>P03995 Related reagents</p> <p>Q28115 Related reagents</p> <p>Entrez Gene:</p> <p>2670 GFAP Related reagents</p> <p>24387 Gfap Related reagents</p> <p>281189 GFAP Related reagents</p> <p>14580 Gfap Related reagents</p>
RRID	AB_10613491
Specificity	<p>Mouse anti Bovine GFAP antibody, clone 4A11 recognizes glial fibrillary acidic protein, or GFAP, a class III intermediate filament. During development of the central nervous system, GFAP distinguishes astrocytes from other glial cells. It is involved in various cellular functions, such as cell structure and movement, cell communication and the functioning of the blood-brain barrier. It also plays a role in mitosis by adjusting the filament network present in the cell.</p> <p>Defects in GFAP are a cause of Alexander disease, a rare disorder of the central nervous system affecting mostly males. It is a progressive leukoencephalopathy resulting in mental and physical retardation, dementia, seizures and early death.</p> <p>Mouse anti Bovine GFAP antibody, clone 4A11 may be used in conjunction with Mouse anti Bovine GFAP antibody, clone 1B4 (MCA4733) and Mouse anti Bovine GFAP antibody, clone 2E1 (MCA4734) for increased sensitivity when used in immunohistology.</p>
Histology Positive	Brain

Control Tissue

References

1. Pegram, C.N. *et al.* (1985) Monoclonal antibodies reactive with epitopes restricted to glial fibrillary acidic proteins of several species. [Neurochem Pathol. 3 \(2\): 119-38.](#)
2. McLendon, R.E. *et al.* (1986) The immunohistochemical application of three anti-GFAP monoclonal antibodies to formalin-fixed, paraffin-embedded, normal and neoplastic brain tissues. [J Neuropathol Exp Neurol. 45 \(6\): 692-703.](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA4684>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M367887:200529'

Printed on 19 Oct 2023