

## Datasheet: MCA4694

<b>Description:</b>	RAT ANTI MOUSE IgA HEAVY CHAIN
<b>Specificity:</b>	IgA HEAVY CHAIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	LO-MA-10
<b>Isotype:</b>	IgM
<b>Quantity:</b>	0.25 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified IgM - liquid
<b>Preparation</b>	Purified IgM from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgM concentration 1.0mg/ml

**External Database****Links****UniProt:**

[P01878](#)   [Related reagents](#)

**Entrez Gene:**

[16061](#)   Igh-VJ558   [Related reagents](#)

**RRID**

AB\_1605055

**Specificity**

**Rat anti Mouse IgA Heavy Chain antibody, clone LO-MA-10** recognises the alpha heavy chain of mouse immunoglobulin. It does not react with human IgG or IgM. It also does not react with chicken, rabbit, goat, sheep, bovine, horse, dog, pig or baboon IgG.

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Guarantee**

12 months from date of despatch

**Health And Safety Information**

Material Safety Datasheet documentation #10040 available at:  
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

**Regulatory**

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rat IgM (302001...) [Biotin](#), [RPE](#)

Goat Anti Rat IgM (STAR116...) [HRP](#)

**North & South**   Tel: +1 800 265 7376

**America**   Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M384226:210513'

**Printed on 21 Mar 2022**

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)