

Datasheet: MCA4699F

BATCH NUMBER 0509R

Description:	MOUSE ANTI HUMAN CD266:FITC
Specificity:	CD266
Other names:	TWEAK-R
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	ITEM-4
Isotype:	IgG2b
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid

Max Ex/Em

Fluorophore	Excitation Max (nm)	Emission Max (nm)
FITC	490	525

Preparation

Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution	Phosphate buffered saline pH7.4
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1mg/ml
Immunogen	Human CD266 transfected P815 cells.
External Database Links	<p>UniProt: Q9NP84 Related reagents</p> <p>Entrez Gene: 51330 TNFRSF12A Related reagents</p>
Synonyms	FN14
RRID	AB_10851349
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the P3U1 mouse myeloma cell line.
Specificity	<p>Mouse anti Human CD266 antibody, clone ITEM-4 recognizes human Tumor necrosis factor receptor superfamily member 12A, also known as CD266, Fibroblast growth factor-inducible immediate-early response protein 14 or Tweak-receptor. CD266 is a 102 amino acid ~14 kDa type I transmembrane glycoprotein and member of the tumour necrosis factor receptor superfamily (TNFRSF12A) expressed by a variety of tissues, including the heart, lung and kidney, and at elevated levels by human brain and liver tumors.</p> <p>The interaction between CD266 and its ligand CD255 (TWEAK), induces the proliferation and migration of endothelial cells and also angiogenesis, and studies have shown that CD255-induced apoptosis of some tumor cell lines is mediated by CD266 (Nakayama <i>et al.</i> 2003).</p> <p>Mouse anti Human CD266 antibody, clone ITEM-4 is cross-reactive with mouse CD266, and also with Rat-1 fibroblastoid cells in flow cytometry.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Nakayama, M. <i>et al.</i> (2003) Fibroblast growth factor-inducible 14 mediates multiple pathways of TWEAK-induced cell death. J Immunol. 170 (1): 341-8. 2. Nakayama, M. <i>et al.</i> (2003) Characterization of murine TWEAK and its receptor (Fn14) by monoclonal antibodies. Biochem Biophys Res Commun. 306 (4): 819-25. 3. Winkles, J.A. and Brown,S. (2010) Targeted Therapy of Fn14-Positive Breast Tumors Using a TWEAK-CytotoxinFusion Protein or Noncovalent Complex. U.S. Army Medical Research and Materiel Command, Annual Report 4. Girgenrath, M. <i>et al.</i> (2006) TWEAK, via its receptor Fn14, is a novel regulator of

mesenchymal progenitor cells and skeletal muscle regeneration. [EMBO J. 25: 5826-39.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10041 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA4699F>
10041

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL:FITC \(MCA691F\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M367911:200529'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)