

Datasheet: MCA4699GA

Description:	MOUSE ANTI HUMAN CD266
Specificity:	CD266
Other names:	TWEAK-R
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	ITEM-4
Isotype:	IgG2b
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/25 - 1/200
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Human CD266 transfected P815 cells.
External Database Links	<p>UniProt: Q9NP84 Related reagents</p> <p>Entrez Gene: 51330 TNFRSF12A Related reagents</p>
Synonyms	FN14
RRID	AB_10891057
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the P3U1 mouse myeloma cell line.
Specificity	<p>Mouse anti Human CD266 antibody, clone ITEM-4 recognizes human tumor necrosis factor receptor superfamily member 12A, also known as CD266, Fibroblast growth factor-inducible immediate-early response protein 14 or Tweak-receptor. CD266 is a 102 amino acid ~14 kDa type I transmembrane glycoprotein and member of the tumour necrosis factor receptor superfamily (TNFRSF12A) expressed by a variety of tissues, including the heart, lung and kidney, and at elevated levels by human brain and liver tumors.</p> <p>The interaction between CD266 and its ligand CD255 (TWEAK), induces the proliferation and migration of endothelial cells and also angiogenesis, and CD255-induced apoptosis of some tumor cell lines is mediated by CD266 (Nakayama et al. 2003).</p> <p>Mouse anti Human CD266 antibody, clone ITEM-4 is cross-reactive with mouse CD266, and also with Rat-1 fibroblastoid cells in flow cytometry.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Nakayama, M. <i>et al.</i> (2003) Fibroblast growth factor-inducible 14 mediates multiple pathways of TWEAK-induced cell death. J Immunol. 170 (1): 341-8. 2. Nakayama, M. <i>et al.</i> (2003) Characterization of murine TWEAK and its receptor (Fn14) by monoclonal antibodies. Biochem Biophys Res Commun. 306 (4): 819-25. 3. Winkles, J.A. and Brown,S. (2010) Targeted Therapy of Fn14-Positive Breast Tumors Using a TWEAK-CytotoxinFusion Protein or Noncovalent Complex. U.S. Army Medical Research and Materiel Command, Annual Report

4. Girgenrath, M. *et al.* (2006) TWEAK, via its receptor Fn14, is a novel regulator of mesenchymal progenitor cells and skeletal muscle regeneration. [EMBO J. 25: 5826-39.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA4699GA>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M384086:210513'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)