

Datasheet: MCA4701

Description:	HAMSTER ANTI MOUSE CD27
Specificity:	CD27
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	LG.3A10
Isotype:	IgG
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/10 - 1/20
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	<p>Reacts with: Human, Rat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 0.5mg/ml
Immunogen	Armenian hamster fibroblast line ARHO12 transfected with mouse CD27 cDNA.
External Database Links	<p>UniProt:</p> <p>P41272 Related reagents</p> <p>P26842 Related reagents</p> <p>Q501W2 Related reagents</p> <p>Entrez Gene:</p> <p>21940 Cd27 Related reagents</p> <p>939 CD27 Related reagents</p>
Synonyms	Tnfrsf7, TNFRSF7
RRID	AB_2073423
Fusion Partners	Armenian hamster spleen cells from immunized animals were fused with cells of the murine SP2/0 non-secreting myeloma.
Specificity	<p>Hamster anti Mouse CD27 antibody, clone LG.3A10 recognizes mouse, rat and human CD27, a member of the tumor necrosis factor receptor superfamily. It is expressed in the thymus and spleen, but not in non-lymphoid tissues. CD27 acts as a ligand for CD70, which is expressed in T and B cells, and acts as a co-stimulatory signal in lymphocyte proliferation and immunoglobulin production.</p> <p>CD27 lacks a death domain in its short cytoplasmic tail, but can induce apoptosis by binding with the apoptosis regulatory protein Siva (also known as CD27-binding protein).</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Gravestain, L.A. <i>et al.</i> (1996) CD27 cooperates with the pre-T cell receptor in the regulation of murine T cell development. J Exp Med. 184 (2): 675-85. 2. Sakanishi, T. & Yagita, H. (2010) Anti-tumor effects of depleting and non-depleting anti-CD27 monoclonal antibodies in immune-competent mice. Biochem Biophys Res Commun. 393 (4): 829-35.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA4701>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Hamster IgG (STAR104...) [DyLight@800](#), [FITC](#)

Goat Anti Hamster IgG (STAR79...) [Biotin](#), [FITC](#), [HRP](#)

Recommended Negative Controls

[HAMSTER \(ARMENIAN\) IgG NEGATIVE CONTROL \(MCA2356\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M390180:210825'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)