

Datasheet: MCA4733GA

BATCH NUMBER 170902

Description:	MOUSE ANTI BOVINE GFAP
Specificity:	GFAP
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	1B4
Isotype:	IgG2b
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			1/200 - 1/1000
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1/200 - 1/2000
Immunofluorescence			▪	
Radioimmunoassays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires protein digestion pre-treatment of paraffin sections e.g. trypsin or pronase.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Human, Mouse, Rat, Sheep, Dog, Pig, Rabbit, Guinea Pig, Chicken
N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Bovine spinal cord homogenate.

External Database Links

- UniProt:**
- [P14136](#) [Related reagents](#)
 - [P47819](#) [Related reagents](#)
 - [P03995](#) [Related reagents](#)
 - [Q28115](#) [Related reagents](#)

- Entrez Gene:**
- [2670](#) GFAP [Related reagents](#)
 - [24387](#) Gfap [Related reagents](#)
 - [281189](#) GFAP [Related reagents](#)
 - [14580](#) Gfap [Related reagents](#)

RRID AB_2021036

Specificity **Mouse anti Bovine GFAP antibody, clone 1B4** recognizes glial fibrillary acidic protein, or GFAP, a class III intermediate filament. During development of the central nervous system, GFAP distinguishes astrocytes from other glial cells. It is involved in various cellular functions, such as cell structure and movement, cell communication and the functioning of the blood-brain barrier. It also plays a role in mitosis by adjusting the filament network present in the cell.

Defects in GFAP are a cause of Alexander disease, a rare disorder of the central nervous system affecting mostly males. It is a progressive leukoencephalopathy resulting in mental and physical retardation, dementia, seizures and early death.

Mouse anti Bovine GFAP antibody, clone 1B4 may be used in conjunction with clones [4A11](#) and [2E1](#) for increased sensitivity when used in immunohistology.

Histology Positive Control Tissue Brain

References

1. Pegram, C.N. *et al.* (1985) Monoclonal antibodies reactive with epitopes restricted to glial fibrillary acidic proteins of several species. [Neurochem Pathol. 3 \(2\): 119-38.](#)
2. McLendon, R.E. *et al.* (1986) The immunohistochemical application of three anti-GFAP monoclonal antibodies to formalin-fixed, paraffin-embedded, normal and neoplastic brain tissues. [J Neuropathol Exp Neurol. 45 \(6\): 692-703.](#)
3. Schwabenlander, M. *et al.* (2016) Brain, Craniofacial, and Dental Lesions of a Free-ranging Gray Wolf (*Canis lupus*) Implicated in a Human Attack in Minnesota, USA. [J Wildl Dis. 52 \(1\): 131-7.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA4733GA>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M438129:250327'

