

## Datasheet: MCA4737

|                      |  |
|----------------------|--|
| <b>Description:</b>  | MOUSE ANTI CLOSTRIDIUM DIFFICILE TOXIN B |
| <b>Specificity:</b>  | CLOSTRIDIUM DIFFICILE TOXIN B            |
| <b>Format:</b>       | Purified                                 |
| <b>Product Type:</b> | Monoclonal Antibody                      |
| <b>Clone:</b>        | 5156                                     |
| <b>Isotype:</b>      | IgG1                                     |
| <b>Quantity:</b>     | 0.1 mg                                   |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                            | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry             |     |    | ▪              |                    |
| Immunohistology - Frozen   |     |    | ▪              |                    |
| Immunohistology - Paraffin |     |    | ▪              |                    |
| ELISA                      | ▪   |    |                | 1/20 - 1/200       |
| Immunoprecipitation        |     |    | ▪              |                    |
| Western Blotting           |     |    | ▪              |                    |
| Immunofluorescence         |     |    | ▪              |                    |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

|                                 |   |
|---------------------------------|---|
| <b>Target Species</b>           | Bacterial   |
| <b>Product Form</b>             | Purified IgG - liquid   |
| <b>Preparation</b>              | Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant |
| <b>Buffer Solution</b>          | Phosphate buffered saline   |
| <b>Preservative Stabilisers</b> | 0.1% Sodium Azide (NaN <sub>3</sub> )   |
| <b>Approx. Protein</b>          | IgG concentration 0.1mg/ml  |

## Concentrations

---

### External Database

#### Links

#### UniProt:

[P18177](#)

[Related reagents](#)

---

#### Synonyms

tcdB

---

#### RRID

AB\_1658139

---

#### Specificity

**Mouse anti Clostridium difficile Toxin B antibody, clone 5156** recognizes *Clostridium difficile* toxin B. No reaction is observed with toxin A. *Clostridium difficile* is a gram-positive motile bacterium which is the leading cause of diarrhoea in developed countries. The incidence of disease in humans varies greatly with age, spore density and the administration of antibiotics.

Toxin B (TcdB) is a large exotoxin. Its role in disease is less well understood than that of toxin A, though it seems clear that TcdB contributes to disease. Both toxins modify the Ras superfamily of small GTPases via glycosylation, inactivating them and leading to the disruption of vital signaling pathways in the cell.

---

#### Further Reading

1. Voth, D.E. & Ballard, J.D. (2005) Clostridium difficile toxins: mechanism of action and role in disease. [Clin Microbiol Rev. 18 \(2\): 247-63.](#)

---

#### Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

#### Guarantee

12 months from date of despatch

---

#### Health And Safety Information

Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA4737>  
10040

---

#### Regulatory

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...) [HRP](#)  
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)  
Goat Anti Mouse IgG (STAR76...) [RPE](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Rabbit Anti Mouse IgG (STAR9...)

[FITC](#)

Goat Anti Mouse IgG (H/L) (STAR117...)

[Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

Goat Anti Mouse IgG (STAR70...)

[FITC](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M421857:230731'

**Printed on 12 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)