

Datasheet: MCA4773T

Description:	MOUSE ANTI HUMAN CHROMOGRANIN A
Specificity:	CHROMOGRANIN A
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	PHE5
Isotype:	IgG1
Quantity:	25 µg

Product Details

RRID AB_10671592

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1/250 - 1/1000
ELISA			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species Human

Species Cross Reactivity Does not react with:Rat

Product Form Purified IgG - liquid

Preparation Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide (NaN₃)

Carrier Free Yes

Approx. Protein IgG concentration 1.0mg/ml

Concentrations

Immunogen Human pheochromocytoma tissue.

External Database Links

UniProt:

[P10645](#) [Related reagents](#)

Entrez Gene:

[1113](#) CHGA [Related reagents](#)

Specificity

Mouse anti Human Chromogranin A antibody, clone PHE5 specifically recognizes chromogranin A (CgA), a secretory glycoprotein and member of the chromogranin/secretogranin family, widely expressed by cells of the neuroendocrine system and related tumours, which is used in immunohistochemistry as a valuable diagnostic marker for pheochromocytomas, paragangliomas and carcinoid, pancreatic and prostate tumours.

CgA constitutes one of the most abundant components of secretory granules, and is the precursor to several functional peptides including vasostatin and pancreastatin.

Mouse anti Human Chromogranin A antibody, clone PHE5 is reported as suitable for use in western blotting, identifying identical bands to clone [LK2H10](#).

Histology Positive Control Tissue

Normal Pancreas

References

1. Bordi, C. *et al.* (1988) Comparative study of seven neuroendocrine markers in pancreatic endocrine tumours. [Virchows Arch A Pathol Anat Histopathol. 413 \(5\): 387-98.](#)

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)

Goat Anti Mouse IgG (STAR77...) [HRP](#)

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Goat Anti Mouse IgG (STAR76...) [RPE](#)

Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@680](#),
[DyLight@800](#), [FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M333997:181127'

Printed on 11 Oct 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)