

Datasheet: MCA490

BATCH NUMBER 150338

Description:	MOUSE ANTI RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Specificity:	RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Other names:	RSV
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	RSV3216 (B016)
Isotype:	IgG2b
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA (1)	▪			
Immunoprecipitation	▪			
Western Blotting	▪			
Immunofluorescence	▪			
Immuno-electron Microscopy	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)Suitable as capture reagent with MCA491 as detection reagent in sandwich ELISA (See [Adams et al. 2010](#) for details).

Target Species	Viral
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on protein A from tissue culture supernatant.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Bovine RSV strains: 127, SNK and 9007. Human RSV strains: Long, Randall, 8/60, and A/2.
RRID	AB_2231368
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
Specificity	<p>Mouse anti respiratory syncytial virus fusion protein antibody, clone RSV3216 recognises an epitope within the RSV fusion protein (46 kDa and 22 kDa s-s linked glycoprotein).</p> <p>Mouse anti respiratory syncytial virus fusion protein antibody, clone RSV3216 can be used in IFA studies in conjunction with MCA491G (clone RSV3132).</p>
Purity	>90% IgG
References	<ol style="list-style-type: none"> Mason, S.W. <i>et al.</i> (2004) Polyadenylation-dependent screening assay for respiratory syncytial virus RNA transcriptase activity and identification of an inhibitor. Nucleic Acids Res. 32 (16): 4758-67. Adams, O. <i>et al.</i> (2010) Palivizumab-resistant human respiratory syncytial virus infection in infancy. Clin Infect Dis. 2010 Jul 15;51(2):185-8. Eckardt-Michel, J. <i>et al.</i> (2008) The fusion protein of respiratory syncytial virus triggers p53-dependent apoptosis. J Virol. 82: 3236-49. Kolokoltsov, A.A. <i>et al.</i> (2007) Small interfering RNA profiling reveals key role of clathrin-mediated endocytosis and early endosome formation for infection by respiratory syncytial virus. J Virol. 81: 7786-800. Riffault, S. <i>et al.</i> (2006) Replication of bovine respiratory syncytial virus in murine cells depends on type I interferon-receptor functionality. J Gen Virol. 87: 2145-8. Schlender, J. <i>et al.</i> (2005) Inhibition of toll-like receptor 7- and 9-mediated alpha/beta interferon production in human plasmacytoid dendritic cells by respiratory syncytial virus and measles virus. J Virol. 79 (9): 5507-15. Schlender, J. <i>et al.</i> (2003) Respiratory syncytial virus (RSV) fusion protein subunit F2, not attachment protein G, determines the specificity of RSV infection. J Virol. 77: 4609-16. Zimmer, G. <i>et al.</i> (2005) A chimeric respiratory syncytial virus fusion protein functionally replaces the F and HN glycoproteins in recombinant Sendai virus. J Virol. 79: 10467-77. Zimmer, G. <i>et al.</i> (2003) Virokinin, a bioactive peptide of the tachykinin family, is released from the fusion protein of bovine respiratory syncytial virus. J Biol Chem. 278: 46854-61. Riffault, S. <i>et al.</i> (2010) A new subunit vaccine based on nucleoprotein nanoparticles

- confers partial clinical and virological protection in calves against bovine respiratory syncytial virus. [Vaccine. 28: 3722-34.](#)
11. Yunus, A. *et al.* (2010) Elevated temperature triggers human respiratory syncytial virus F protein six-helix bundle formation. [Virology. 396: 226-37.](#)
 12. Zimmer, G. *et al.* (2002) Cleavage at the furin consensus sequence RAR/KR(109) and presence of the intervening peptide of the respiratory syncytial virus fusion protein are dispensable for virus replication in cell culture. [J Virol. 76 \(18\): 9218-24.](#)
 13. Adams, O. *et al.* (2013) Genetic analysis and antigenic characterization of human respiratory syncytial virus group A viruses isolated in Germany 1996-2008. [Virus Genes. 47 \(2\): 210-8.](#)
 14. Zimmer, G. *et al.* (2001) Proteolytic activation of respiratory syncytial virus fusion protein. Cleavage at two furin consensus sequences. [J Biol Chem. 276 \(34\): 31642-50.](#)
 15. Zimmer, G. *et al.* (2001) N-glycans of F protein differentially affect fusion activity of human respiratory syncytial virus. [J Virol. 75 \(10\): 4744-51.](#)
 16. Shaw, C.A. *et al.* (2013) The role of non-viral antigens in the cotton rat model of respiratory syncytial virus vaccine-enhanced disease. [Vaccine. 31 \(2\): 306-12.](#)
 17. Sandhu, J.S. *et al.* (2000) Oral immunization of mice with transgenic tomato fruit expressing respiratory syncytial virus-F protein induces a systemic immune response. [Transgenic Res. 9 \(2\): 127-35.](#)
 18. Hägglund S *et al.* (2014) Characterization of an experimental vaccine for bovine respiratory syncytial virus. [Clin Vaccine Immunol. 21 \(7\): 997-1004.](#)
 19. Hervé, P.L. *et al.* (2016) RSV N-nanorings fused to palivizumab-targeted neutralizing epitope as a nanoparticle RSV vaccine. [Nanomedicine. Aug 20. pii: S1549-9634\(16\)30114-9. \[Epub ahead of print\]](#)
 20. Fleming, E.H. *et al.* (2006) Respiratory syncytial virus F envelope protein associates with lipid rafts without a requirement for other virus proteins. [J Virol. 80 \(24\): 12160-70.](#)
 21. Wong TM *et al.* (2014) Respiratory syncytial virus (RSV) infection in elderly mice results in altered antiviral gene expression and enhanced pathology. [PLoS One. 9 \(2\): e88764.](#)
 22. Bird, G.H. *et al.* (2014) Mucosal delivery of a double-stapled RSV peptide prevents nasopulmonary infection. [J Clin Invest. 124 \(5\): 2113-24.](#)
 23. Sampayo-Escobar, V. *et al.* (2018) Osteopontin plays a pivotal role in increasing severity of respiratory syncytial virus infection. [PLoS One. 13 \(4\): e0192709.](#)

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA490</p> <p>10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M368063:200529'

Printed on 12 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)