

Datasheet: MCA493

BATCH NUMBER 158671

Description:	MOUSE ANTI hCG (ALPHA 6 EPITOPE)
Specificity:	hCG (ALPHA 6 EPITOPE)
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	INN-hCG-72
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			
Radioimmunoassays	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein	Isotype concentration 4.0 mg/ml

Concentrations

Immunogen Human CG.

External Database Links

UniProt:

[P01215](#) [Related reagents](#)

Entrez Gene:

[1081](#) CGA [Related reagents](#)

RRID AB_321601

Fusion Partners Spleen cells from immunised mice were fused with cells of the X63-Ag8.653 mouse myeloma cell line.

Specificity **Mouse anti Human hCG antibody, clone INN-hCG-72** recognises exclusively free hCG alpha. It reacts with the alpha 6 epitope on hCG alpha. No cross reaction is observed with holo-hCG in direct binding RIA of ELISA.

References

1. Berger, P. *et al.* (1990) Monoclonal antibodies against the free subunits of human chorionic gonadotrophin. [J Endocrinol. 125 \(2\): 301-9.](#)
 2. Madersbacher, S. *et al.* (1992) Free alpha-subunit, free beta-subunit of human chorionic gonadotropin (hCG), and intact hCG in sera of healthy individuals and testicular cancer patients. [Clin Chem. 38 \(3\): 370-6.](#)
 3. Berger, P. *et al.* (1993) Variants of human chorionic gonadotropin from pregnant women and tumor patients recognized by monoclonal antibodies. [J Clin Endocrinol Metab. 77 \(2\): 347-51.](#)
 4. Madersbacher, S. *et al.* (1993) Serum glycoprotein hormones and their free alpha-subunit in a healthy elderly population selected according to the SENIEUR protocol. Analyses with ultrasensitive time resolved fluoroimmunoassays. [Mech Ageing Dev. 71 \(3\): 223-33.](#)
 5. Dirnhofer, S. *et al.* (1994) Free alpha subunit of human chorionic gonadotropin: molecular basis of immunologically and biologically active domains. [J Endocrinol. 140 \(1\): 145-54.](#)
 6. Hiro'oka, T. *et al.* (2000) Disulfide bond mutations in follicle-stimulating hormone result in uncoupling of biological activity from intracellular behavior. [Endocrinology. 141 \(12\): 4751-6.](#)
-

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA493>
10040

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M384241:210513'

Printed on 18 Jan 2024