

Datasheet: MCA54R

BATCH NUMBER 170836

Description:	MOUSE ANTI RAT CD43
Specificity:	CD43
Other names:	LEUKOSIALIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	W3/13
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/25 - 1/100
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.

Target Species	Rat
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Rat thymocyte membrane glycoproteins.
External Database Links	<p>UniProt: P13838 Related reagents</p> <p>Entrez Gene: 24796 Spn Related reagents</p>
RRID	AB_2194201
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
Specificity	<p>Mouse anti Rat CD43 antibody, clone W3/13 recognizes the rat CD43 cell surface antigen, also known as leukosialin, sialophorin or W3/13 antigen. CD43 is a 371 amino acid ~95 kDa heavily glycosylated single pass type 1 transmembrane glycoprotein (Killeen et al. 1987) expressed by all leucocytes with the exception of B lymphocytes. CD43, in mice acts as a T-cell counter-receptor for CD169 (Siglec-1) suggesting a role in cell-cell interactions (van den Berg et al. 2001)</p> <p>Mouse anti Rat CD43 antibody, clone W3/13 is routinely tested in flow cytometry on rat splenocytes.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Brown, W.R.A. <i>et al.</i> (1981) Identification of a glycoprotein-like molecule at the cell surface of rat thymocytes. Nature. 289: 456-460. 2. Barclay, A. N. (1981) The localization of populations of lymphocytes defined by monoclonal antibodies in rat lymphoid tissues. Immunology. 42: 593-600 3. Forbes, R. D. C. <i>et al.</i> (1983) Leukocyte subsets in first-set rat cardiac allograft rejection. A serial immunohistologic study using monoclonal antibodies Transplantation. 36: 681-686 4. Jung, S. <i>et al.</i> (1994) Therapeutic effect of transforming growth factor-beta 2 on actively induced EAN but not adoptive transfer EAN Immunology. 83: 545-551. 5. Bataller, R. <i>et al.</i> (2003) Prolonged infusion of angiotensin II into normal rats induces stellate cell activation and proinflammatory events in liver. Am J Physiol Gastrointest Liver Physiol. 285: G642-651 6. Rice, E.K. <i>et al.</i> (2003) Induction of MIF synthesis and secretion by tubular epithelial

- cells: a novel action of angiotensin II. [Kidney Int. 63 \(4\): 1265-75.](#)
7. Schwab, J.M. *et al.* (2005) Spinal cord injury induces early and persistent lesional P2X4 receptor expression. [J Neuroimmunol. 163: 185-9.](#)
 8. Conrad, S. *et al.* (2005) Prolonged lesional expression of RhoA and RhoB following spinal cord injury. [J Comp Neurol. 487: 166-75.](#)
 9. Zhang, Z. *et al.* (2008) FTY720 ameliorates experimental autoimmune neuritis by inhibition of lymphocyte and monocyte infiltration into peripheral nerves. [Exp Neurol. 210: 681-90.](#)
 10. Yamanaka, Y. *et al.* (2011) Immunohistochemical analysis of subcutaneous tissue reactions to methacrylate resin-based root canal sealers. [Int Endod J. 44: 669-75.](#)
 11. Dort, J. *et al.* (2012) Beneficial effects of cod protein on skeletal muscle repair following injury. [Appl Physiol Nutr Metab. 37 \(3\): 489-98.](#)
 12. Duchesne, E. *et al.* (2013) Mast cells can regulate skeletal muscle cell proliferation by multiple mechanisms. [Muscle Nerve. 48 \(3\): 403-14.](#)
 13. Xu, K. *et al.* (2016) Expression of aryl hydrocarbon receptor in rat brain lesions following traumatic brain injury. [Diagn Pathol. 11 \(1\): 72.](#)
 14. Dort, J. *et al.* (2016) Shrimp Protein Hydrolysate Modulates the Timing of Proinflammatory Macrophages in Bupivacaine-Injured Skeletal Muscles in Rats. [Biomed Res Int. 2016: 5214561.](#)
 15. Zhang, Z.M. *et al.* (2016) Lesional accumulation of CD8(+) cells in sciatic nerves of experimental autoimmune neuritis rats. [Neurol Sci. 37 \(2\): 199-203.](#)
 16. Ornellas, F.M. *et al.* (2019) Mesenchymal Stromal Cells Induce Podocyte Protection in the Puromycin Injury Model. [Sci Rep. 9 \(1\): 19604.](#)
 17. Grad, E. *et al.* (2019) Monocyte Modulation by Liposomal Alendronate Improves Cardiac Healing in a Rat Model of Myocardial Infarction [Regenerative Engineering and Translational Medicine. 5 \(3\): 280-9.](#)
 18. Kami, K. *et al.* (2019) Role of 72-kDa Heat Shock Protein in Heat-stimulated Regeneration of Injured Muscle in Rat. [J Histochem Cytochem. 67 \(11\): 791-9.](#)
 19. Cakała-Jakimowicz, M. & Puzianowska-Kuznicka, M. (2022) Towards Understanding the Lymph Node Response to Skin Infection with Saprophytic *Staphylococcus epidermidis*. [Biomedicines. 10 \(5\): 1021.](#)
 20. Merlini, A. *et al.* (2022) Distinct roles of the meningeal layers in CNS autoimmunity. [Nat Neurosci. 25 \(7\): 887-99.](#)

Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA54R
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA1209\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M438123:250326'

Printed on 31 May 2026