

Datasheet: MCA554XZ

**BATCH NUMBER 153806**

<b>Description:</b>	MOUSE ANTI HUMAN CD35:Preservative Free
<b>Specificity:</b>	CD35
<b>Other names:</b>	CR1
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	E11
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1)This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.Sodium citrate buffer pH 6.0 is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Cynomolgus monkey, Rhesus Monkey, Baboon <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1mg/ml
<b>Immunogen</b>	Human acute monocytic leukaemia cells.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P17927</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">1378</a>    CR1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	C3BR
<b>RRID</b>	AB_323143
<b>Fusion Partners</b>	Spleen cells from immunized BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
<b>Specificity</b>	<p><b>Mouse anti Human CD35 antibody, clone E11</b> recognizes complement receptor type 1, also known as CD35 or the C3b/C4b receptor. CD35 a single chain, type I single pass transmembrane glycoprotein which exists in four allotypic forms (A,B,C,D) of ~190 kDa, ~220 kDa, ~160 kDa and ~250 kDa respectively (<a href="#">Khera et al. 2009</a>).</p> <p>CD35 is expressed by granulocytes, monocytes, B cells and some T cells. Functionally, CD35 acts as a receptor for the C3b component of complement.</p>
<b>Flow Cytometry</b>	Use 10ul of the diluted antibody to label 10 <sup>6</sup> cells or 100ul whole blood.
<b>Histology Positive Control Tissue</b>	Tonsil
<b>References</b>	<ol style="list-style-type: none"> <li>Birmingham, D.J. <i>et al.</i> (1996) The baboon erythrocyte complement receptor is a glycoposphatidylinositol-linked protein encoded by a homologue of the human CR1-like genetic element. <a href="#">J Immunol. 157 (6): 2586-92.</a></li> <li>Sopper, S. <i>et al.</i> (1997) Lymphocyte subsets and expression of differentiation markers in blood and lymphoid organs of rhesus monkeys. <a href="#">Cytometry. 29 (4): 351-62.</a></li> </ol>

3. Weissgerber P *et al.* (2003) Investigation of mechanisms involved in phagocytosis of *Legionella pneumophila* by human cells. [FEMS Microbiol Lett. 219 \(2\): 173-9.](#)
4. Schmidt, C.Q. *et al.* (2016) Selectivity of C3-opsonin targeted complement inhibitors: A distinct advantage in the protection of erythrocytes from paroxysmal nocturnal hemoglobinuria patients. [Immunobiology. 221 \(4\): 503-11.](#)
5. Noda, G.S. *et al.* (2020) Specificities and isotypes of erythrocytes autoantibodies in patients with warm autoimmune hemolytic anemia [Rev Cubana Hematol Inmunol Hemoter 36\(4\): e1283.](#)
6. Klaas, M. *et al.* (2023) Sialoadhesin (CD169/Siglec-1) is an extended molecule that escapes inhibitory cis-interactions and synergizes with other macrophage receptors to promote phagocytosis. [Glycoconj J. 40 \(2\): 213-23.](#)
7. Wagstaffe, H.R. *et al.* (2023) The Development of Immunological Assays to Evaluate the Level and Function of Antibodies Induced by *Klebsiella pneumoniae* O-Antigen Vaccines. [mSphere. : e0068022.](#)

---

**Storage**

Store at -20°C only.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10162 available at: <https://www.bio-rad-antibodies.com/SDS/MCA554XZ>  
10162

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

- |   |   |
|---|---|
| Rabbit Anti Mouse IgG (STAR12...)       | <a href="#">RPE</a>   |
| Goat Anti Mouse IgG IgA IgM (STAR87...) | <a href="#">HRP</a>   |
| Goat Anti Mouse IgG (STAR76...)         | <a href="#">RPE</a>   |
| Goat Anti Mouse IgG (STAR70...)         | <a href="#">FITC</a>  |
| Goat Anti Mouse IgG (H/L) (STAR117...)  | <a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> ,<br><a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> ,<br><a href="#">FITC</a> , <a href="#">HRP</a> |
| Goat Anti Mouse IgG (STAR77...)         | <a href="#">HRP</a>   |
| Goat Anti Mouse IgG (Fc) (STAR120...)   | <a href="#">FITC</a> , <a href="#">HRP</a>  |
| Rabbit Anti Mouse IgG (STAR13...)       | <a href="#">HRP</a>   |
| Rabbit Anti Mouse IgG (STAR9...)        | <a href="#">FITC</a>  |

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M368236:200529'

**Printed on 25 Apr 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)