

Datasheet: MCA5641P

BATCH NUMBER 160254

Description:	MOUSE ANTI LIPOPROTEIN LIPASE:HRP
Specificity:	LIPOPROTEIN LIPASE
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	5D2
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoprecipitation	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Species Cross Reactivity	<p>Reacts with: Rat, Human, Chicken, Guinea Pig</p> <p>Does not react with: Mouse</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Purified bovine milk lipoprotein lipase.
External Database Links	<p>UniProt: P11151 Related reagents</p> <p>Entrez Gene: 280843 LPL Related reagents</p>
RRID	AB_10963722
Specificity	<p>Mouse anti lipoprotein lipase antibody, clone 5D2 recognizes an epitope within amino acids 380-410 of lipoprotein lipase (LPL), a member of the AB hydrolase superfamily, which plays a pivotal role in lipoprotein metabolism and transport, acting as the key enzyme in the hydrolysis of triglycerides and very low density lipoproteins (VLDLs), and the release of free fatty acids into peripheral tissues.</p> <p>A deficiency of LPL can result in hypertriglyceridemia, and LPL plays a critical role in the pathogenesis of atherosclerosis, and in particular the relationship between LPL and apolipoprotein E (ApoE), both of which are secreted in significant amounts by macrophages in developing arterial wall lesions.</p> <p>Mouse anti Lipoprotein lipase antibody, clone 5D2 is a unique antibody which differentiates between monomeric inactive and dimeric active LPL, and binds to LPL sequences involved in LPL, LPL receptor, and heparin interactions. Mouse anti Lipoprotein lipase antibody, clone 5D2 inhibits the activity of human LPL (Chang <i>et al.</i> 1998).</p>
ELISA	MCA5641P can be used in a direct ELISA, or as the detection reagent in a sandwich ELISA with MCA5641 or MCA5641GA as the capture antibody.
Western Blotting	MCA5641P detects a band of approximately 53kDa using partially purified LPL from postheparin plasma.
References	<ol style="list-style-type: none"> Peterson, J. <i>et al.</i> (1992) Human lipoprotein lipase: relationship of activity, heparin affinity, and conformation as studied with monoclonal antibodies. J Lipid Res. 33 (8): 1165-70. Chang, S.F. <i>et al.</i> (1998) Detailed characterization of the binding site of the lipoprotein lipase-specific monoclonal antibody 5D2. J Lipid Res. 39 (12): 2350-9. Hussain, M.M. <i>et al.</i> (2000) High affinity binding between lipoprotein lipase and lipoproteins involves multiple ionic and hydrophobic interactions, does not require enzyme activity, and is modulated by glycosaminoglycans. J Biol Chem. 275 (38): 29324-30.

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10094 available at: <https://www.bio-rad-antibodies.com/SDS/MCA5641P>
10094

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M393577:220112'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)