

## Datasheet: MCA5669

|                      |                       |
|----------------------|-----------------------|
| <b>Description:</b>  | MOUSE ANTI HUMAN CD15 |
| <b>Specificity:</b>  | CD15                  |
| <b>Other names:</b>  | LEWIS X               |
| <b>Format:</b>       | Purified              |
| <b>Product Type:</b> | Monoclonal Antibody   |
| <b>Clone:</b>        | MCS-1                 |
| <b>Isotype:</b>      | IgG3                  |
| <b>Quantity:</b>     | 0.2 mg                |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                            | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry             | ▪   |    |                |                    |
| Immunohistology - Frozen   |     |    | ▪              |                    |
| Immunohistology - Paraffin |     |    | ▪              |                    |
| ELISA                      |     |    | ▪              |                    |
| Immunoprecipitation        |     |    | ▪              |                    |
| Western Blotting           |     |    | ▪              |                    |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

|                                       |   |
|---------------------------------------|---|
| <b>Target Species</b>                 | Human   |
| <b>Product Form</b>                   | Purified IgG - liquid   |
| <b>Buffer Solution</b>                | Phosphate buffered saline   |
| <b>Preservative Stabilisers</b>       | <0.1% Sodium Azide (NaN <sub>3</sub> )<br>1% Bovine Serum Albumin |
| <b>Approx. Protein Concentrations</b> | IgG concentration 0.1mg/ml  |

RRID AB\_10672021

**Specificity** **Mouse anti Human CD15 antibody, clone MCS-1** recognizes human CD15, otherwise known as Lewis X and stage-specific embryonic antigen-1 (SSEA-1). The CD15 antigen is itself a terminal trisaccharide, Gal-beta1-4 (Fuc-alpha1-3) GlcNAc (3-FAL), which is found on a variety of glycoproteins and glycolipids. CD15 is predominately expressed by peripheral blood granulocytes and various epithelial tissues. It plays an important role in cell adhesion and has been used as a marker for diagnosis of Hodgkin's lymphoma.

**Flow Cytometry** Use 20ul of the suggested working dilution to label  $1 \times 10^6$  cells in 100ul.

**References** 1. Oxford University Press. (1998) . Leucocyte Typing IV. p798-810.

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Guarantee** 12 months from date of despatch

**Health And Safety Information** Material Safety Datasheet documentation #10041 available at: 10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

**Regulatory** For research purposes only

## Related Products

### Recommended Secondary Antibodies

|   |   |
|---|---|
| Goat Anti Mouse IgG (STAR77...)         | <a href="#">HRP</a>   |
| Rabbit Anti Mouse IgG (STAR12...)       | <a href="#">RPE</a>   |
| Human Anti Mouse IgG3 (HCA039...)       | <a href="#">FITC</a> , <a href="#">HRP</a>  |
| Goat Anti Mouse IgG (STAR70...)         | <a href="#">FITC</a>  |
| Goat Anti Mouse IgG IgA IgM (STAR87...) | <a href="#">Alk. Phos.</a> , <a href="#">HRP</a>  |
| Rabbit Anti Mouse IgG (STAR9...)        | <a href="#">FITC</a>  |
| Goat Anti Mouse IgG (STAR76...)         | <a href="#">RPE</a>   |
| Rabbit Anti Mouse IgG (STAR13...)       | <a href="#">HRP</a>   |
| Goat Anti Mouse IgG (Fc) (STAR120...)   | <a href="#">FITC</a> , <a href="#">HRP</a>  |
| Goat Anti Mouse IgG (H/L) (STAR117...)  | <a href="#">Alk. Phos.</a> , <a href="#">DyLight@488</a> , <a href="#">DyLight@550</a> ,<br><a href="#">DyLight@650</a> , <a href="#">DyLight@680</a> , <a href="#">DyLight@800</a> ,<br><a href="#">FITC</a> , <a href="#">HRP</a> |

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

**Printed on 14 Sep 2022**

---

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)