

# Datasheet: MCA5750Z

**BATCH NUMBER 163905**

<b>Description:</b>	MOUSE ANTI HUMAN VEGF:Preservative Free
<b>Specificity:</b>	VEGF
<b>Other names:</b>	VPF
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	4B2-8
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			2.0 - 4.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.25 - 0.50ug/ml
Functional Assays	▪			0.1 - 0.2ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - lyophilized
<b>Reconstitution</b>	Reconstitute with 0.5 ml distilled water
	Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture

supernatant

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml after reconstitution
<b>Immunogen</b>	<a href="#">Highly pure recombinant human VEGF<sub>165</sub></a> (PHP091).
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P15692</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">7422</a> VEGFA <a href="#">Related reagents</a>
<b>Synonyms</b>	VEGF
<b>RRID</b>	AB_10897188
<b>Specificity</b>	<p><b>Mouse anti Human Vascular endothelial growth factor antibody, clone 4B2-8</b> recognizes human Vascular endothelial growth factor A , also known as VEGF-A or Vascular permeability factor. VEGF-A is a 206 amino acid ~45 kDa disulphide linked homodimeric heparin-binding cytokine and member of the PDGF/VEGF growth factor family, which exists as several splice variants of the VEGF gene, and is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis.</p> <p>The functional forms of VEGF are classified according to the number of amino acid residues within a constitutive polypeptide chain e.g VEGF<sub>165</sub>, and all VEGF family members act through binding to the tyrosine kinase receptors VEGFR-1 (Flt-1) and VEGFR-2 (KDR/CD309), whilst the closely-related VEGF-C and VEGF-D, can also bind to VEGFR-3 (Flt-4).</p>
<b>ELISA</b>	MCA5750Z may be used in a sandwich ELISA with <a href="#">AHP831B</a> as the detection antibody and <a href="#">PHP091</a> as the standard.
<b>Western Blotting</b>	MCA5750Z detects a band of approximately 38.2kDa using recombinant human VEGF <sub>165</sub> <a href="#">PHP091</a> under non-reducing conditions.
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.</p> <p>After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not</p>

recommended. Avoid repeated freezing and thawing as this may denature the antibody.  
Should this product contain a precipitate we recommend microcentrifugation before use.

<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10294 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA5750Z">https://www.bio-rad-antibodies.com/SDS/MCA5750Z</a> 10294
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">Alk. Phos.</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M402158:220718'

Printed on 18 Jan 2024