

Datasheet: MCA5750Z

BATCH NUMBER 163905

Description:	MOUSE ANTI HUMAN VEGF:Preservative Free
Specificity:	VEGF
Other names:	VPF
Format:	Preservative Free
Product Type:	Monoclonal Antibody
Clone:	4B2-8
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			2.0 - 4.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.25 - 0.50ug/ml
Functional Assays	▪			0.1 - 0.2ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.5 ml distilled water
	Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture

supernatant

Buffer Solution Phosphate buffered saline

Preservative Stabilisers None present

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0 mg/ml after reconstitution

Immunogen [Highly pure recombinant human VEGF₁₆₅](#) (PHP091).

External Database Links

UniProt:

[P15692](#)

[Related reagents](#)

Entrez Gene:

[7422](#)

VEGFA

[Related reagents](#)

Synonyms VEGF

RRID AB_10897188

Specificity

Mouse anti Human Vascular endothelial growth factor antibody, clone 4B2-8

recognizes human Vascular endothelial growth factor A , also known as VEGF-A or Vascular permeability factor. VEGF-A is a 206 amino acid ~45 kDa disulphide linked homodimeric heparin-binding cytokine and member of the PDGF/VEGF growth factor family, which exists as several splice variants of the VEGF gene, and is a potent stimulator of angiogenesis of both normal and cancerous cells, and acts as a regulator of vasculogenesis.

The functional forms of VEGF are classified according to the number of amino acid residues within a constitutive polypeptide chain e.g VEGF₁₆₅, and all VEGF family members act through binding to the tyrosine kinase receptors VEGFR-1 (Flt-1) and VEGFR-2 (KDR/CD309), whilst the closely-related VEGF-C and VEGF-D, can also bind to VEGFR-3 (Flt-4).

ELISA

MCA5750Z may be used in a sandwich ELISA with [AHP831B](#) as the detection antibody and [PHP091](#) as the standard.

Western Blotting

MCA5750Z detects a band of approximately 38.2kDa using recombinant human VEGF₁₆₅ [PHP091](#) under non-reducing conditions.

Storage

Prior to reconstitution store at -20°C.

After reconstitution store at -20°C.

This product should be stored undiluted. Storage in frost-free freezers is not

recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10294 available at: https://www.bio-rad-antibodies.com/SDS/MCA5750Z 10294
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M402158:220718'

Printed on 12 Aug 2023