

Datasheet: MCA5790F

Description:	RAT ANTI MOUSE TIM-3:FITC
Specificity:	TIM-3
Other names:	T-CELL IMMUNOGLOBULIN MUCIN 3
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	RMT3-23
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse						
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid						
Max Ex/Em	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>FITC</td> <td>490</td> <td>525</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	FITC	490	525
Fluorophore	Excitation Max (nm)	Emission Max (nm)					
FITC	490	525					
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant						
Buffer Solution	Phosphate buffered saline						
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin						
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml						

Immunogen	Tim-3-Ig, consisting of amino acids 1-191 of mouse Tim-3, coupled with the Fc portion of mouse IgG2a.
External Database Links	<p>UniProt: Q8VIM0 Related reagents</p> <p>Entrez Gene: 171285 Havcr2 Related reagents</p>
Synonyms	Tim3, Timd3
Fusion Partners	Lymph node cells from immunised SD rats were fused with cells of the P3U1 myeloma cell line.
Specificity	<p>Rat anti Mouse Tim-3 antibody, clone RMT3-23 recognizes mouse T-cell immunoglobulin mucin 3, a transmembrane glycoprotein and member of the immunoglobulin superfamily, preferentially expressed by differentiated CD4+ Th1 cells and CD8+ Tc1 (cytotoxic) cells, but not by CD4+ Th2 cells or CD8+ Tc2 cells.</p> <p>Tim-3 acts as a negative regulator of Th1-mediated inflammatory diseases including GVHD, type I diabetes, and EAE (experimental autoimmune encephalomyelitis), promotes immunological tolerance, and is a receptor for galectin-9, shown to induce cell death of Th1 cells <i>in vitro</i>. Tim-3 is also involved in the regulation of macrophage activation, acts as a phagocytic receptor, and has been shown to recognize apoptotic cells via the IgV domain.</p> <p>Rat anti Mouse Tim-3 antibody, clone RMT3-23 significantly inhibits the efficient phagocytosis of apoptotic cells by peritoneal exudate Mac1+ cells (Nakayama et al. 2009).</p>
Flow Cytometry	<p>Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.</p> <p>The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors. This may be reduced by using SeroBlock FcR (BUF041A or BUF041B).</p>
References	<ol style="list-style-type: none"> Oikawa, T. <i>et al.</i> (2006) Preferential involvement of Tim-3 in the regulation of hepatic CD8+ T cells in murine acute graft-versus-host disease. J Immunol. 177 (7): 4281-7. Nakayama, M. <i>et al.</i> (2009) Tim-3 mediates phagocytosis of apoptotic cells and cross-presentation. Blood. 113 (16): 3821-30.
Further Reading	1. Freeman, G.J. <i>et al.</i> (2010) TIM genes: a family of cell surface phosphatidylserine receptors that regulate innate and adaptive immunity. Immunol Rev. 235 (1): 172-89.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. This product is photosensitive and should be</p>

protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:FITC \(MCA1212F\)](#)

Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)

[MOUSE SEROBLOCK FcR \(BUF041B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M368406:200529'

Printed on 22 Mar 2021