

Datasheet: MCA5862G

BATCH NUMBER 0112R

Description:	MOUSE ANTI HUMAN ALPHA FETOPROTEIN
Specificity:	ALPHA FETOPROTEIN
Other names:	AFP
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	125.6G7
Isotype:	IgG1
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			5µg/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)

Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
External Database Links	<p>UniProt: P02771 Related reagents</p> <p>Entrez Gene: 174 AFP Related reagents</p>
Synonyms	HPAFP
RRID	AB_10672456
Specificity	<p>Mouse anti Human alpha fetoprotein, clone 125.6G7 is a monoclonal antibody recognizing alpha fetoprotein (AFP), a member of the albuminoid protein superfamily, along with albumin and vitamin D binding protein. Members of the albuminoid superfamily comprise molecules with three characteristic globular domains (Mizejewski 1995) .</p> <p>AFP is a ~70 kDa plasma protein found in fetuses over four weeks old, produced by the yolk sac and the fetal liver, with the highest concentration reached during the 12th to 16th week of gestation. AFP acts as a transport protein for various ligands such as bilirubin, fatty acids, steroids and heavy metals. After birth, plasma AFP levels drop rapidly with only trace amounts remaining in healthy adults. Changes in plasma levels of AFP have been documented in a variety of cancers including hepatocellular cancer (Sato <i>et al.</i> 1993)</p> <p>Clone 125.6G7 shows <0.1% reactivity with human serum albumin and IgG. Strong expression of AFP has been noted in the hepatic cell line HepG2.</p>
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA5862G 10040</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M368453:200529'

Printed on 15 Jul 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)