

Datasheet: MCA5932P

BATCH NUMBER 1012

Description:	MOUSE ANTI CHICKEN INTERLEUKIN-6:HRP
Specificity:	IL-6
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	AV129
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/20,000 - 1/200,000
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chicken
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.01% Thiomersal
Stabilisers	HRP Stabiliser (BUF052A)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen Chicken interleukin-6 (IL-6)

External Database

Links

UniProt:

[Q90YI0](#) [Related reagents](#)

Entrez Gene:

[395337](#) IL6 [Related reagents](#)

Fusion Partners

Spleen cells from immunised mice were fused with cells of the SP2/0 myeloma cell line

Specificity

Mouse anti Chicken Interleukin-6 antibody, clone AV129 recognizes chicken interleukin-6 (IL-6), a cytokine that can function in either a pro-inflammatory or in an anti-inflammatory capacity.

Chicken interleukin-6 is produced by a range of cell types including macrophages, endothelial cells, T-cells and fibroblasts. IL-6 is a multipotential cytokine with the ability to induce proliferation and differentiation in B-cells and T-cells and acts as an important mediator in acute phase response.

This antibody was raised as part of the UK BBSRC funded avian immunology toolbox project and is one of a number reagents now available from Bio-Rad developed from this project.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10131 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA5932P>
10131

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[RECOMBINANT CHICKEN INTERLEUKIN-6 \(PAP003\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M368487:200529'

Printed on 12 Aug 2023

