

Datasheet: MCA5933GA

BATCH NUMBER 171506

Description:	MOUSE ANTI SHEEP CD335
Specificity:	CD335
Other names:	Gr13.1, NKp46
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	EC1.1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/25 - 1/200
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Mouse anti Sheep CD335 antibody, clone EC1.1 is suitable for use on paraffin sections when fixed with zinc salts Results may be enhanced by the use of a using TSA (Tyramide Signal Amplification) biotin signal enhancement system.

Target Species	Sheep
Species Cross Reactivity	<p>Reacts with: Bovine, Goat</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>

Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	OvNKp46-mFcy2b fusion protein
External Database Links	<p>UniProt:</p> <p>E5LDC2 Related reagents</p> <p>Q863H2 Related reagents</p> <p>Entrez Gene:</p> <p>369024 NCR1 Related reagents</p>
Fusion Partners	Spleen cells from immunized Balb/c mice were fused with cells of the X63 myeloma cell line
Specificity	<p>Mouse anti Sheep CD335, clone EC1.1, recognizes the ovine homologue of human CD335, also known as NKp46 and natural cytotoxicity triggering receptor 1. CD335 is a type I transmembrane protein, with two extracellular C2-type immunoglobulin-like domains, which function as an activating receptor.</p> <p>In mammals, the CD335 antigen is confined to NK cells (Westergard et. al. 2004) and is primarily expressed by a subset of natural killer (NK) cells which are CD3^{-ve}, CD4^{-ve} and CD2^{+ve} with a small population also being CD8^{+ve} (Storset et. al. 2004). CD335 is not detected on B cells, T cells, monocytes or granulocytes. NK cells are cytotoxic and form a key component of the innate immune system involved in mediating immunity to tumors and pathogens in particular the process where the host recognizes and rejects tumor and virally infected cells.</p> <p>Mouse anti Sheep CD335, clone EC1.1 (originally named as Gr13.1) is the first antibody specifically raised to ovine CD335 and provides a reagent enabling identification of ovine NKp46 positive NK cells (Connelley et. al. 2011). In addition to clone EC1.1, clones AKS1, which recognizes CD335 (NKp46) in bovine and other ruminants, and VIV-KM1, which recognizes porcine CD335, are also available from Bio-Rad.</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul

References

- Olsen, L. *et al.* (2013) Characterization of NCR1+ cells residing in lymphoid tissues in the gut of lambs indicates that the majority are NK cells [Vet Res. 44:109.](#)
- Goh, S. *et al.* (2016) Identification of *Theileria lestoquardi* Antigens Recognized by CD8+ T Cells. [PLoS One. 11 \(9\): e0162571.](#)
- Kolar, Q.K. *et al.* (2020) Anatomical distribution of respiratory tract leukocyte cell subsets in neonatal calves. [Vet Immunol Immunopathol. 227: 110090.](#)
- Marzo, S. *et al.* (2021) Characterisation of dendritic cell frequency and phenotype in bovine afferent lymph reveals kinetic changes in costimulatory molecule expression [Vet Immunol Immunopathol. 19 Nov: 110363.](#)

Further Reading

- Storset, A. K. *et al.* (2004) NKp46 defines a subset of bovine leukocytes with natural killer cell characteristics. [Eur J Immunol. 34:669-676](#)
- Westgaard, I.H. *et al.* (2004) Rat NKp46 activates natural killer cell cytotoxicity and is associated with FcεpsilonRIγ and CD3zeta. [J Leukoc Biol. 76 \(6\): 1200-6.](#)
- Elh mouzi-Younes, J. *et al.* (2010) Ovine CD16+/CD14- blood lymphocytes present all the major characteristics of natural killer cells. [Vet Res. 41:4](#)
- Boysen, P. *et al.* (2011) Natural killer cells in free-living *Mus musculus* have a primed phenotype. [Mol Ecol. 20 \(23\): 5103-10.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements The Renaissance TSA Biotin system is a registered trademark of PerkinElmer, Waltham, USA

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA5933GA>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

- Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
- Goat Anti Mouse IgG (STAR70...) [FITC](#)
- Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
- Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
- Goat Anti Mouse IgG (STAR77...) [HRP](#)
- Goat Anti Mouse IgG (STAR76...) [RPE](#)
- Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Recommended Useful Reagents

[MOUSE ANTI HUMAN CD16 \(MCA5665GA\)](#)
[MOUSE ANTI HUMAN CD16:FITC \(MCA5665F\)](#)
[MOUSE ANTI HUMAN CD16:RPE \(MCA5665PE\)](#)
[MOUSE ANTI SHEEP CD16-Like \(MCA919GA\)](#)
[RAT ANTI HUMAN CD3:Pacific Blue® \(MCA1477PB\)](#)
[RAT ANTI HUMAN CD3:FITC \(MCA1477F\)](#)
[MOUSE ANTI BOVINE CD205 \(MCA1651GA\)](#)
[MOUSE ANTI BOVINE CD205:FITC \(MCA1651F\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M404488:220830'

Printed on 11 Nov 2025

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)