

Datasheet: MCA6071F

Description:	MOUSE ANTI HUMAN CD193:FITC
Specificity:	CD193
Other names:	CCR3
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	5E8
Isotype:	IgG2b
Quantity:	100 TESTS/0.5ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide (NaN ₃)		
Stabilisers	0.2% Bovine Serum Albumin		

External Database Links

UniProt:

[P51677](#)

[Related reagents](#)

Entrez Gene:

[1232](#) CCR3 [Related reagents](#)

Synonyms

CMKBR3

Specificity

Mouse anti Human CD193 antibody, clone 5E8 recognizes the human CD193 cell surface antigen, also known as CC-chemokine receptor 3 (CCR3).

CD193 is a G protein-coupled receptor which binds a number of chemokines including eotaxin, RANTES and MCP-1. It is highly expressed on basophils, eosinophils and airway epithelial cells and is known to be involved in the development of allergic diseases. CD193 is also known to be an entry co-receptor for HIV-1.

Flow CytometryUse 5ul of the suggested working dilution to label 1×10^6 cells in 100ul**Storage**

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL:FITC \(MCA691F\)](#)

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M344218:190111'

Printed on 29 Aug 2021