

## Datasheet: MCA6089

**BATCH NUMBER 151959**

<b>Description:</b>	HUMAN ANTI VEGF (BEVACIZUMAB BIOSIMILAR)
<b>Specificity:</b>	VEGF
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	Bevacizumab
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.2 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Human IgG1 (kappa) antibody - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.02% Proclin 300
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml
<b>Immunogen</b>	Human VEGF

**Specificity** **Human Anti-VEGF antibody** is a research grade biosimilar of the monoclonal antibody drug bevacizumab. It is a recombinant human IgG1 kappa antibody with variable regions

from the therapeutic bevacizumab antibody clone.

This antibody is for research use only. It is suitable as a reference standard in a pharmacokinetic (PK) bridging ELISA with anti-bevacizumab antibodies, for example [HCA182](#) and [HCA185](#). It can also be used as capture and detection reagent in an anti-drug antibody (ADA) bridging ELISA, with an anti-bevacizumab antibody in IgG1 format as a reference standard, for example [HCA184](#) or [HCA185](#).

Bevacizumab, also known as Avastin, is a humanized monoclonal antibody (IgG1 kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A stimulates new blood vessel formation in certain diseases, especially in cancer. Avastin is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

This product is NOT FOR THERAPEUTIC USE.

[View a summary of supporting anti-bevacizumab antibody products.](#)

---

<b>Purity</b>	>98% by SDS PAGE
---------------	------------------

---

<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
----------------	--

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10199 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA6089">https://www.bio-rad-antibodies.com/SDS/MCA6089</a> 10199
--------------------------------------	---

---

<b>Licensed Use</b>	For in vitro research purposes only. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
---------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA182\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA184\)](#)

[HUMAN ANTI BEVACIZUMAB:HRP \(HCA184P\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA185\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M357840:190905'

**Printed on 01 Mar 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)