

Datasheet: MCA6089

BATCH NUMBER 153789

Description:	HUMAN ANTI VEGF (BEVACIZUMAB BIOSIMILAR)
Specificity:	VEGF
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	Bevacizumab
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Human IgG1 (kappa) antibody - liquid
Preparation	Purified recombinant IgG prepared by affinity chromatography on Protein A from a mammalian cell line
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Proclin 300
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Human VEGF
Specificity	Human Anti-VEGF antibody is a research grade biosimilar of the monoclonal antibody

drug bevacizumab. It is a recombinant human IgG1 kappa antibody with variable regions from the therapeutic bevacizumab antibody clone.

This antibody is for research use only. It is suitable as a reference standard in a pharmacokinetic (PK) bridging ELISA with anti-bevacizumab antibodies, for example [HCA182](#) and [HCA185](#). It can also be used as capture and detection reagent in an anti-drug antibody (ADA) bridging ELISA, with an anti-bevacizumab antibody in IgG1 format as a reference standard, for example [HCA184](#) or [HCA185](#).

Bevacizumab, also known as Avastin, is a humanized monoclonal antibody (IgG1 kappa) that inhibits angiogenesis by binding to VEGF-A. VEGF-A stimulates new blood vessel formation in certain diseases, especially in cancer. Avastin is used to treat various cancers, including colorectal, lung, breast, glioblastoma, kidney and ovarian.

This product is NOT FOR THERAPEUTIC USE.

[View a summary of supporting anti-bevacizumab antibody products.](#)

Purity	>98% by SDS PAGE
---------------	------------------

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #20463 available at: https://www.bio-rad-antibodies.com/SDS/MCA6089 20463
--------------------------------------	--

Licensed Use	For in vitro research purposes only. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
---------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK001P\)](#)

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA182\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA184\)](#)

[HUMAN ANTI BEVACIZUMAB:HRP \(HCA184P\)](#)

[HUMAN ANTI BEVACIZUMAB \(HCA185\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M373141:200827'

Printed on 01 Mar 2024
