

Datasheet: MCA6099F

BATCH NUMBER 172564

Description:	MOUSE ANTI PIG B CELL SUBSET:FITC
Specificity:	B CELL SUBSET
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	BB6-10A10
Isotype:	IgM
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig		
Product Form	Purified IgM conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgM prepared by combination of precipitation and chromatography techniques		
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)		
Approx. Protein Concentrations	IgM concentration 0.5 mg/ml		
Specificity	Mouse anti Pig B cell subset antibody, clone BB6-10A10, appears to recognize an		

immature pig B cell population. Clone BB6-10A10 stains the following percentages of surface immunoglobulin (slg) positive cells in porcine lymphoid tissues: 11-15% in peripheral blood, 22-23% in mesenteric lymph nodes, 79-87% in ileal Peyer's patches, 6-17 in spleen, and 0% in thymus ([Boersma et al. 2001](#)).

Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul
References	1. Boersma, W.J. et al. (2001) Summary of workshop findings for porcine B-cell markers. Vet Immunol Immunopathol. 80 (1-2): 63-78.
Storage	This product is shipped at ambient temperature. Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA6099F
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE ANTI PIG CD45RA:RPE \(MCA1751PE\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M441107:250523'

Printed on 29 Jan 2026