

Datasheet: MCA6099GA

BATCH NUMBER 164898

Description:	MOUSE ANTI PIG B CELL SUBSET
Specificity:	B CELL SUBSET
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	BB6-10A10
Isotype:	IgM
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified IgM - liquid
Preparation	Purified IgM prepared by combination of precipitation and chromatography techniques
Buffer Solution	Borate buffered saline
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgM concentration 0.5 mg/ml

Specificity **Mouse anti Pig B cell subset antibody, clone BB6-10A10**, appears to recognize an immature pig B cell population. Clone BB6-10A10 stains the following percentages of surface immunoglobulin (sIg) positive cells in porcine lymphoid tissues: 11-15% in peripheral blood, 22-23% in mesenteric lymph nodes, 79-87% in ileal Peyer's patches,

6-17 in spleen, and 0% in thymus ([Boersma et al. 2001](#)).

Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul
-----------------------	------------------------------------------------------------------------------------

References	1. Boersma, W.J. <i>et al.</i> (2001) Summary of workshop findings for porcine B-cell markers. Vet Immunol Immunopathol. 80 (1-2): 63-78.
-------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10077 available at: https://www.bio-rad-antibodies.com/SDS/MCA6099GA 10077
--------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgM (STAR138...) [Alk. Phos.](#)
Human Anti Mouse IgM (HCA040...) [FITC](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgM NEGATIVE CONTROL \(MCA692\)](#)

Recommended Useful Reagents

[MOUSE ANTI PIG CD45RA \(MCA1751GA\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M392060:211020'

Printed on 19 Jan 2024