

Datasheet: MCA6101

BATCH NUMBER 167582

Description:	HUMAN ANTI CD52 (ALEMTUZUMAB BIOSIMILAR)
Specificity:	CD52
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	Campath-1H
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Human IgG1 (kappa) antibody - liquid
Preparation	Purified recombinant IgG prepared by affinity chromatography on Protein A from a mammalian cell line
Source	HEK
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% ProClin 300
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Human lymphocytes

External Database Links	UniProt: P31358 Related reagents
	Entrez Gene: 1043 CD52 Related reagents
Synonyms	CDW52, HE5
Specificity	<p>Human Anti-CD52 Antibody, clone Campath-1H is a research grade biosimilar of the monoclonal antibody drug alemtuzumab. It is a recombinant human IgG1 kappa antibody with variable regions from the therapeutic antibody Campath-1H.</p> <p>Human Anti-CD52 Antibody is for research use only. This product is supplied in a research grade format, i.e. in PBS buffer with preservative, it is not manufactured or formulated in the same way as the therapeutic reference product. A preservative free version is available on request.</p> <p>It is suitable as a reference standard in a pharmacokinetic (PK) bridging ELISA with anti-alemtuzumab antibodies, for example HCA198 and HCA175. It can also be used as capture and detection reagent in an anti-drug antibody (ADA) bridging ELISA, with an anti-alemtuzumab antibody in IgG1 format as reference standard, for example HCA175.</p> <p>Alemtuzumab is a drug used in the treatment of chronic lymphocytic leukemia (CLL), cutaneous T-cell lymphoma (CTCL) and T-cell lymphoma. Alemtuzumab is also used in some conditioning regimes for bone marrow transplantation, kidney transplantation and islet cell transplantation along with clinical trial protocols for treatment of certain autoimmune diseases including multiple sclerosis.</p> <p>This product is NOT FOR THERAPEUTIC USE.</p> <p>View a summary of supporting anti-alemtuzumab antibody products</p>
Purity	>98% by SDS PAGE
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	12 months from date of despatch
Acknowledgements	ProClin is a trademark of The Dow Chemical Company ("Dow") or an affiliated company of Dow.
Health And Safety	Material Safety Datasheet documentation #20463 available at:

Information <https://www.bio-rad-antibodies.com/SDS/MCA610120463>

Licensed Use For in vitro research purposes only. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[LYNX RAPID HRP ANTIBODY CONJUGATION KIT \(LNK002P\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M423211:231006'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)